

GARD™skin Assay Protocol

Genomic Allergen Rapid Detection (GARD) assay for assessment of skin sensitizers (GARDskin) assay protocol

The Genomic Allergen Rapid Detection skin (GARDskin) is an *in vitro* assay that provides binary hazard identification of skin sensitizers (i.e. UN GHS Category 1) versus non-sensitizers). The method evaluates the transcriptional patterns of a specific genomic biomarker signature, the GARDskin prediction signature (GPS), in the SenzaCell cell line exposed to test chemicals.

Résumé

The purpose of the test method is to contribute to the identification of skin sensitizers and non-sensitizers by providing information on the chemical-induced alteration of gene expression of mechanistically relevant genomic biomarkers associated with the activation of dendritic cells (DC). The activation process through which DC change from antigen processing to antigen presenting cells addresses the third key event of the skin sensitisation Adverse Outcome Pathway (AOP).

Experimental Description

Endpoint Measurement

After test chemical exposure of SenzaCell cells, the expression of the GPS, comprising 196 gene transcripts, is quantified utilizing the nanoString technology.

Endpoint Value(s)

The gene expression of the GPS is analysed with a fixed SVM algorithm that generates "Decision Values" (DVs). From triplicate replicate experiments the mean DV is calculated. If the mean DV is ≥ 0 the test chemical is classified as a skin sensitizer. If the mean DV is < 0 the test chemical is classified as a non-sensitizer.

Experimental System

The biological Test System is a human myeloid dendritic like cell line, SenzaCell, which is a subclone from the MUTZ-3 cell line. SenzaCells are deposited at ATCC and can be purchased from SenzaGen AB.

Discussion

GARDskin assay overview




Figure 1 GARDskin assay workflow. First the SenzaCells are cultured according to standard procedures. Then the Test Item is solubilized and the Test Item cytotoxicity is assessed. Main stimulations follow to generate RNA samples which are the basis for the endpoint measurement, data analysis and classification.

Figure 1 describes the workflow of the GARDskin assay, starting with routine cell culturing of SenzaCell cell line. In the GARDskin assay, cell stimulations with a test chemical are performed for two reasons: 1) cytotoxicity assessment to identify a relevant stimulation concentration, and 2) for the *GARDskin Main Stimulation* to harvest RNA.

In the cytotoxicity assessment, the test chemical is screened for cytotoxic effects, to identify an appropriate concentration (i.e. the concentration that yields a Relative viability of ~90%, Rv90) to be used as the input concentration in the *GARDskin Main Stimulation*. The cells are exposed to a range of concentrations of the test chemical, originating from a serial dilution.

Once the input concentration of the test chemical is found, cells are exposed to the test chemical again in the *GARDskin Main Stimulation* with the concentration identified during the cytotoxicity assessment. Main stimulations are repeated three times to achieve three biological replicate samples. Thus, every GARDskin assessment of a test chemical is based on three replicate *GARDskin Main Stimulation*s.

The endpoint measurements of GARDskin, i.e. the quantification of the GARDskin biomarker signature mRNA transcripts, is performed on total RNA purified from cells from the *GARDskin Main Stimulation*. The quantification is performed using the NanoString nCounter instrument and endpoint specific biomarker CodeSets. The result is analyzed with the *GARD Data Analysis Application* (GDAA) and the test chemical is classified as a sensitizer or non-sensitizer.

Status

Development

Lund University (Sweden)
SenzaGen AB (Sweden)

Known Laboratory Use

EuroFins GB
BRT
MB Research laboratory
SenzaGen AB

Participation in Validation Studies

SenzaGen has organized and completed an inter-laboratory validation study in which the method was transferred to 2 laboratories (see below) and the predictive capacity was also evaluated in each laboratory (Johansson et al. Toxicological Sciences, 2019).

Burleson Research Technologies Inc, 120 First Flight Lane, Morrisville, NC, USA.a.
Eurofins BioPharma Product Testing Munich GmbH, 82152 Planegg/Munich Germany

An independent peer review by the EURL ECVAM Scientific Advisory Committee (ESAC) regarding the SenzaGen-coordinated Performance Standards-based validation of the GARDskin method for skin sensitisation testing has been published in July 2021 (ESAC, 2021).

Regulatory Acceptance

The test method is under review by the OECD (2021).

Proprietary and/or Confidentiality Issues

Intellectual property rights

The SenzaCell cell line is available under a license agreement upon request. The IP rights of the GARD biomarker signatures and any assay utilizing the signatures, in its entirety or parts thereof, are owned by SenzaGen AB.

This Document, its contents and any other intellectual property rights referred to in this document, are the property of SenzaGen AB and are intended solely for guidance concerning performing the GARDskin assay and shall belong to and vest in SenzaGen AB absolutely to the fullest extent permitted by law.

Any patent, registered trademarks, copyright and/or other intellectual property right of SenzaGen AB shall belong to and vest in SenzaGen AB absolutely to the fullest extent permitted by the applicable laws. To SenzaGen knowledge, none of the offered products sold, nor any process or know-how used, by SenzaGen infringes or is alleged to infringe any patent or other intellectual property rights of any other Person.

This document and translations of it may be copied and furnished to others, and derivative works that comment on or otherwise explain it or assist in its implementation may be prepared, copied, published and distributed, in whole or in part, without restriction of any kind, provided that the copyright notice "©2019 SenzaGen AB. All rights reserved " and this paragraph are included on all such copies and derivative works. However, this document itself may not be modified in any way, such as by removing the copyright notice.

Disclaimer

To the maximum extent permitted by law, SenzaGen (on behalf of itself and SenzaGen affiliates, licensors, distributors, agents, consultants, supplier, subcontractors and employees) expressly disclaim all representations and warranties concerning the scope or validity of the intellectual property rights, and expressly disclaim any warranty that the contents of any documents or the developed products referred to herein in connection with any add on content or the use of the intellectual property rights by licensee will not infringe upon patent, copyright, trademark or other proprietary rights of a third party. Any warranty that may be provided in any applicable laws, statute or regulation governing commercial activities is explicitly disclaimed.

Health and Safety Issues

The human myeloid leukemia cell line (SenzaCell) is a cell line of Biosafety level I. As such, no extraordinary safety issues are considered necessary, beyond those considered common for sterile work with mammalian cell lines in laboratories dedicated for such purposes.

The testing laboratory should have in place routines for risk assessment of chemicals. This routine should take into account the information in the SDS and state safe ways of working and waste handling of chemical substances.

The positive control, proficiency chemicals, and vehicles used in the GARDskin assay need to be handled according to available information, such as SDS, since these can be e.g. sensitizing, acutely toxic, and environmentally hazardous etc. The safe handling of a test chemical in the GARDskin assay need to be assessed on a case to case basis. In the case of an unknown test chemical it is recommended to consider it as a sensitizer and a highly toxic compound and use appropriate protection.

The TRIzol or equivalent reagents are corrosive and carcinogenic and should be handled according to instructions provided by the supplier and according to available safety data sheet (SDS).

Propidium iodide is suspected to cause genetically defects and is a known irritant compound which should be used with care according to SDS.

General precautions

Always wear protective clothing and gloves and work in a fume hood when handling chemical substances and the TRIzol reagent. Wear protective glasses and breathing mask when handling the original stocks (powder or liquid) of the chemical substances, and preferably also disposable arm cuffs to avoid contact with the chemical substances.

MSDS Information

Positive control and solvents (negative controls), see CAS number and Catalog number in Table 3, read available MSDS.

Abbreviations and Definitions

Annotation file	A text file used by the digital analyzer during the endpoint analysis to map individual RCC files to substance names and gene signatures.
BSA	Bovine Serum Albumin (used to prepare flow cytometry buffer)
CDF	Cartridge Definition File (defines sample specific data to associate with the data output and defines the parameters for the imaging instrument to use during image collection and processing)
Cell batch	Within the context of this protocol, a unique cell batch is defined as: <ul style="list-style-type: none">- cells originating from different frozen vials, or...- cells originating from the same frozen vial, which have been cultivated separately. A division of cell cultures for the purpose of achieving separate cell batches should be done no sooner than passage 3 after thawing, and no later than at least 2 passages prior to exposure experiments.
DMF	Dimethylformamide (may be used to solubilize test chemicals)
DMSO	Dimethyl Sulfoxide (may be used to solubilize test chemicals)
D-PBS	Dulbecco's Phosphate Buffered Saline (used to prepare flow cytometry wash buffer)
DV	Decision Value (the assay readout)
FBS	Fetal Bovine Serum (a supplement to the cell culture medium)
FITC	Fluorescein isothiocyanate (a derivate of fluorescein which can be conjugated to antibodies. Its fluorescence is detected by flow cytometry analysis)
FSC-A	The Forward Scatter (Area) of the particle population is a parameter that can be detected by a flow cytometry. A higher value corresponds to larger particles.
GARD	Genomic Allergen Rapid Detection
GDAA	GARD Data Analysis Application (a web interfacere used to analyse the GARDskin end-point)
GHS	Globally Harmonized System
GM-CSF	Granulocyte Macrophage Colony Stimulating Factor (a growth factor supplemented to the cell culture medium)
GPS	GARD Prediction Signature (the transcripts measured in the GARDskin assay)

In-well concentration	The test and/or control chemical concentration in the exposed SenzaCell cell culture during a stimulation
mAbs	Monoclonal antibodies (used for the phenotype control)
Main Stimulation	Cell stimulation performed to generate RNA. Each GARDskin test consists of three Main Stimulations
PI	Propidium Iodide (PI is a fluorescent substance that intercalates in DNA. Since viable cells have an intact cell membrane, those will not be stained. The fluorescent signal can be detected by flow cytometry and is used to assess cell viability)
PE	Phycoerythrin (a fluorochrome that can be conjugated to antibodies. Its fluorescence can be detected by flow cytometry analysis).
PPD	p-Phenylenediamine (the GARDskin positive control chemical)
RCC	Reporter Code Count (the file type that is generated by the nCounter System)
RLF	Reporter Library File (a file that contains the code key information used during image processing in the digital analyzer to assign target identities to the barcode)
Rv90	90% relative viability (the concentration of a test chemical or positive control chemical inducing 90 % relative viability. Calculated after flow cytometry analysis of PI stained cell cultures)
SVM	Support Vector Machine
SSC-A	The Side Scatter (Area) of the particle population is a parameter that can be detected by flow cytometry. A higher value corresponds to more granular particles.
TC flask	Tissue Culture flask

PROCEDURE DETAILS

Materials and Preparations

Cell or Experimental system

The human myeloid leukemia cell line SenzaCell is provided by SenzaGen AB and sent to the licensed CRO on dry ice. The provided vial should be stored in liquid nitrogen or a cryogenic freezer. The vial should be expanded and frozen in liquid nitrogen or a cryogenic freezer as a homogenous cell bank according to instructions from the developing laboratory.

Equipment and Consumables

Listed in Table 1-3 are the equipment, consumables, and reagents needed to perform the GARDskin assay.

Table 1 Equipment. Listed are the equipment used at the laboratory of the test method developers. Equipment for which a specific manufacturer is not listed, the source of the equipment is considered arbitrary.

Equipment	Manufacturer
Sterile (LAF) hood for cell culture work, Class II	-
Fume hood for handling of chemicals	-
Heating block for small tubes (operating at 37°C±5°C)	-
Laboratory grade scale, capacity of weighing a minimum of 10 mg with reproducibility	-
Vials for long-term storage of cells in liquid nitrogen	-
2 separate CO ₂ incubators (one for cell line culturing and one for chemical stimulation)	-
Benchtop centrifuge, swing-out rotor, 2-8°C, including adapters for 5/15/50 ml tubes and adapters for 96 well plates	-
Freezer (operating at -18°C to -22°C)	-
Ultra-low freezer (operating at -70°C to -90°C)	-
Refrigerator (operating at 2-8°C)	-
Flow cytometer (minimum equipped with a blue laser, e.g. FACSVerse)*	E.g. BD
Microcentrifuge for 1.5 ml micro tubes	-
Minicentrifuge for 0.2 ml tubes	-
BioAnalyzer 2100*	E.g. Agilent
nCounter MAX, FLEX or sprint system**	NanoString
Thermocycler	-
Centrifuge adapters for 96 well plates	-
Pipetting controller	-
Pipettes 0.1-1000 µl	-
Good to have: Electronic multi-dispenser pipette 1 µl-50 ml	
Good to have: Vacuum system for aspiration of cell supernatant and other liquid waste.	-

* The flow cytometer and BioAnalyzer could be exchanged for an **equivalent instrument**.

**The endpoint analysis is performed with an nCounter instrument, either at the Test Facility or at a NanoString Facility if not available in-house.

Table 2. Consumables Listed are the consumables used at the laboratory of the test method developers. The exchange of any of these articles for an **equivalent product** should not interfere with the protocols and/or results, but needs to be assessed to ensure equivalence, especially the cell culture plastics.

Product	Company*	Catalog Number
TC Flask, 175 cm ²	Corning	431080
TC Flask, 75 cm ²	Corning	430641U
TC Flask, 25 cm ²	Corning	430372
Centrifuge tube, 15 ml	Corning Sarstedt	430791 62.554.502
Centrifuge tube, 50 ml	Corning Sarstedt	430829 62.547.254
12-well plate	Corning	3512
24-well plate	Corning	3524
Stripettes, 5 ml	Corning Sarstedt	4051 86.1253.001
Stripettes, 10 ml	Corning Sarstedt	4101 86.1254.001
Stripettes, 25 ml	Corning Sarstedt	4251 86.1685.001
Stripettes, 50 ml	Corning Sarstedt	4501 86.1256.001
Positive displacement pipette tips, 1 ml	Eppendorf	0030089642
Positive displacement pipette tips, 5 ml	Eppendorf	0030089669
Positive displacement pipette tips, 10 ml	Eppendorf	0030089677
Positive displacement pipette tips, 25 ml	Eppendorf	0030089685
Positive displacement pipette tips, 50 ml	Eppendorf	0030089693
Sterile and RNase-free filter tips 0.1-1000 µl	Sartorius (SafetySpace) Biotix (uTIP)	Several
Cryogenic Vial	Corning	430488
Sample tubes for Flow cytometry	Corning	352052
(Deep 96 well plate if used for flow cytometry)	Corning	3960
0.2 µm sterile filter	Sarstedt	83.1826.001
Syringe Luer, 50 ml	Henke-Sass Wolf	8300006680
1.5 ml micro tubes	Sarstedt	72.690.001
RNase-free 1.5 ml micro tubes	Axygen	311-09-051
(RNase-free 2 ml vials in 96 racks if used for homogenized cell samples)	Micronic	MP42150
RNase-free 0.2 ml tubes	Sarstedt	72.991.002
Nitrile gloves, thickness 0.14 mm	Shieldskin	67625

Table 3. Reagents. The listed producers are made available as a guidance for quality of listed reagents. The exchange of any of these articles for an equivalent product should not interfere with the GARDskin assay, but needs to be assessed to ensure equivalence. The NanoString products, are specifically required in this protocol and equivalent reagents do not, at this time, exist.

Product	Company	Catalog Number
<i>Cell Medium</i>		
MEM/Alpha Modification with L-glut, Ribo-& Deoxyribo	Cytiva	SH30265.01
Fetal Bovine Serum (FBS)** (Gibco)	Life Technologies	10270106
rhGM-CSF (Premium grade. Purity >97%, endotoxin level <0.1 EU/μg cytokine, and activity of ≥5x10 ⁶ IU/mg)	Miltenyi Biotec	130-093-868
<i>Buffers & Solvents</i>		
D-PBS, HyClone or D-PBS, Gibco	Cytiva Life Technologies	SH30028.02 14190144
Bovine Serum Albumin (BSA), Cohn fraction V	Saveen & Werner AB	B2000
TRIzol Reagent or TRI Reagent	Ambion Thermo Fisher Scientific	15596018 AM9738
Ethanol, 95-100%, Udenatured	Solveco	1065
<i>Antibodies & Staining</i>		
Mouse anti-human CD86-FITC	BD	555657
Mouse anti-human HLA-DR-FITC	BD	347400
Mouse anti-human CD34-FITC	BD	555821
Mouse anti-human CD1a-FITC	Agilent Dako	F714101-2
Mouse anti-human CD54-PE	BD	555511
Mouse anti-human CD14-PE	Agilent Dako	R086401-2
Mouse anti-human CD80-PE	BD	340294
Mouse polyclonal anti-IgG1-FITC	BD	555748
Mouse polyclonal anti-IgG1-PE	BD	555749
Propidium Iodide, 50 μg/ml	BD	556463
Trypan Blue Solution, 0.4%	Thermo Scientific	15250061
<i>Reagents & Kits</i>		
Direct-zol RNA MiniPrep	Zymo Research	R2052
RNA 6000 Nano Kit	Agilent	5067-1511
nCounter MAX or Sprint consumables	NanoString	See nanostring product information.
GPS200_v2 CodeSet (GARDskin)	NanoString	(Contact SenzaGen)
<i>Chemicals and solvents</i>		
	<i>CAS no</i>	<i>Catalog no**</i>
p-Phenylenediamine (PPD)	106-50-3	965106
DMSO ≥99.5%	67-68-5	D5879
Acetone	67-64-1	34850
Ethanol, 95-100%	64-17-5	-
Dimethylformamide	68-12-2	1.03053
Isopropanol	67-63-0	I9516
Glycerol	56-81-5	G5516

* Each lot of FBS needs to be assessed for cell culture compatibility.

** Catalog numbers at Sigma-Aldrich for guidance of e.g. purity of each chemical.

Media and Endpoint Assay Solutions

Serum, GM-CSF and antibodies

Prior to performing a GARDskin assay, FBS needs to be assessed and antibodies need to be titrated (see Annex 1. FBS assessment and Annex 2. Antibody titration). The FBS can be aliquoted and stored long-term at -18°C or below. The GM-CSF working stock (150 µg/ml) should be prepared, aliquoted, and stored long-term at -18°C or below. The GM-CSF working stock can be stored short-term at 2-8°C for maximum 1 week.

Cell medium Cell culture medium is prepared in two steps and referred to as follows:

1) Semi-complete medium

MEM/Alpha supplemented with FBS (20%).

Use within 30 days from supplementing FBS. Store at 2-8°C.

2) Complete medium

Semi-complete medium supplemented with GM-CSF (40 ng/ml).

0.26 µl GM-CSF (150 µg/ml) is added per ml semi-complete medium.

Complete medium cannot be stored and should be used directly.

See Table 3 for GM-CSF purity, endotoxin level and activity.

Medium for freezing cells

The SenzaCell cell line is frozen and stored in liquid nitrogen (or in a cryogenic freezer) in complete medium supplemented with DMSO to a final DMSO concentration of 10%.

Flow cytometry Wash buffer

For all washing and staining steps and for suspension of cells prior to flow cytometry run, use flow cytometry wash buffer: D-PBS supplemented with ~0.5-1% (w/w) BSA. The BSA is dissolved in a smaller volume of D-PBS, which is filtered using a 0.2 µm filter prior addition to the D-PBS. The prepared Wash buffer can be stored at 2-8°C for 30 days.

Summary of storage conditions for cell reagents.

Table 4: Storage conditions of routine cell culture and flow cytometry reagents.

Reagent	Storage long term	Storage short term	Assessment/titration
FBS	-18 °C or below until expiry date.	2-8 °C for at least 30 days.	yes
GM-CSF	-18 °C or below until expiry date.	2-8 °C for at least 7 days.	no
Antibodies	2-8 °C until expiry date.	Not applicable	yes
MEM/alpha	2-8 °C until expiry date.	Not applicable	no
Semi-complete medium	2-8 °C for 30 days	2-8 °C for 30 days	no
Complete medium	Not applicable	Use same day as mixing.	no
Flow cytometry wash buffer	Not applicable	2-8 °C for 30 days	no

Flow cytometer instrument setup

Prior to performing the GARDskin assay for the first time, a fluorescence compensation of the flow cytometer should be performed. The compensation should be performed according to the specific instrument. In flow cytometry, "Compensation" is a mathematical correction of a signal overlap between the channels of the emission spectra of different fluorochromes. In the GARDskin assay, the compensation should preferably be performed using the SenzaCell cell line single stained with the mAbs HLA-DR-FITC and CD54-PE. For flow cytometry analysis, appropriate flow rate should be set (FACSVerse 60-120 µl/min).

Controls

For each GARDskin test, relevant controls are analyzed in each of the three replicate *GARD Main Stimulations*. The controls are listed in Table 5, with relevant information for the laboratory work. The GARDskin input concentration of the positive control, p-Phenylenediamine (PPD), should be determined in a cytotoxicity assessment while establishing the method in the laboratory. The input concentration for PPD should be the same independent of operator at the same laboratory.

Table 5. List of controls used in the GARDskin assay.

Substance ID	Control	GARDskin classification	Solvent	GARD Input conc (µM)	Viability
pos ctrl	p-Phenylenediamine (PPD)	Sensitizer	DMSO	As determined in cytotoxicity assessment.	84.5 - 95.4% (Relative Viability)
neg ctrl	Test chemical vehicle*	Non-sensitizer	Not applicable	Concentration corresponding to vehicle concentration in test chemical stimulation well*	≥95.5% (Relative Viability)
unstim ctrl	No stimulant added beside complete cell culture medium	Not applicable	Not applicable	Not applicable	≥84.5% (Absolute Viability)

*See Table 8 on page 18 for available vehicles and corresponding maximum concentration.

Unstimulated controls

The unstimulated control is used for determination of absolute and relative cell viability of cell batches and for normalization purposes in the Data analysis workflow.

Negative controls

The negative control is a vehicle control to verify that cells have not become activated in any steps of the method's experimental procedures, i.e. the vehicle shall not induce cytotoxicity and shall be classified as a non-sensitizer.

Positive controls

The positive control shall verify that the cells are responsive and can become activated upon exposure of a skin sensitizer, i.e. the positive control shall induce cell cytotoxicity and be classified as a sensitizer.

Demonstration of Proficiency

According to OECD TG 442E GARDskin assay laboratories should demonstrate technical proficiency in using the test method prior to routine use of GARDskin. Proficiency is demonstrated by testing of a specified set of proficiency chemicals with known sensitising properties, as listed in Table 6.

Table 6. Substances for demonstrating technical proficiency with GARDskin.

Proficiency substances	CASRN	Reference value ¹	GARDskin prediction
4-nitrobenzyl bromide	100-11-8	Sensitiser	Sensitiser
Propyl gallate	121-79-9	Sensitiser	Sensitiser
Isoeugenol	97-54-1	Sensitiser	Sensitiser
3-(Dimethylamino)-1-propylamine	109-55-7	Sensitiser	Sensitiser
Eugenol	97-53-0	Sensitiser	Sensitiser
Ethylene glycol dimethacrylate	97-90-5	Sensitiser	Sensitiser
Glycerol	56-81-5	Non-sensitiser	Non-Sensitiser
Hexane	110-54-3	Non-sensitiser	Non-Sensitiser
1-Butanol	71-36-3	Non-sensitiser	Non-Sensitiser

¹ NICEATM LLNA database 2010.

Method

Routine Procedures

All cell work should be performed under sterile conditions free of antibiotics; work in a laboratory designed for growth of mammalian cells, use LAF-workbenches and sterile consumables.

All cell centrifugation steps are performed at 300-315xg, 5 min, 2-8°C. All cell incubations are performed in cell incubators at 37°C±1°C and 5%±0.5% CO₂ at saturated humidity.

Cell cultures should not be grown for more than 16 passages (~ 2 months) after thawing. A cell passage is defined for SenzaCell cells as each time the cell culture is counted and split, independently of how the cells has grown i.e. its doubling time (see below sections *Thawing of cells* for details about cell passage numbering, and *Cell seeding for test chemical stimulation* for details about the range of cell passages used in cell stimulation).

For cell maintenance, grow cells in cell culture flasks. For volumes up to 10 ml, use TC Flask 25 cm². For volumes of 10-45 ml, use TC Flask 75 cm². For volumes of 40-120 ml, use TC Flask 162-175 cm². Note that for large cultures, more than one TC 162-175 cm² may be required.

Thawing cells

The SenzaCells are stored in liquid nitrogen (liquid phase) or in a cryogenic freezer, 7 million cells /ml complete medium supplemented with 10% v/v DMSO.

- Thaw the cells by submerging the bottom half of the frozen vial in a ~37°C water.

- Add 10 ml semi-complete medium to a 15 ml tube and transfer the thawed cells to the tube. Centrifuge the cells.
- Remove supernatant by decantation. Resuspend the cell pellet in 5 ml semi-complete medium. Add 0.26 μl GM-CSF per 1 ml of cell suspension to the cell culture.
- Move the cells to a small cell culture flask (TC Flask 25 cm^2) and incubate the cell culture (i.e. cell passage number P0).
- The next day, transfer the cell culture from the cell culture flask to a 50 ml tube. Centrifuge.
- Remove supernatant by decantation. Resuspend in 1 ml semi-complete medium.
- Count the cells. Resuspend the cell culture in semi-complete medium to a volume corresponding to a cell concentration of 2×10^5 cells /ml.
- Add 0.26 μl GM-CSF per 1 ml of cell suspension to the cell culture. Incubate the cell culture (i.e. cell passage number P1).

Cell batches

For convenience, it is preferred to run three cell batches in parallel. Within the GARDassay, a unique cell batch is defined as follows.

- cells originating from different frozen vials, or,
- cells originating from the same frozen vial, which have been cultivated separately. A division of cell cultures for the purpose of achieving separate cell batches should be done no sooner than passage 3 after thawing, and no later than at least 2 passages prior to exposure experiments.

Cell Propagation

Every 3-4 days the cells are counted and propagated to 2×10^5 cells /ml in fresh medium. The cell propagation is preferably performed on Mondays and Thursdays to coincide with cell stimulations (see *Cell seeding for test chemical stimulation*).

- To split the cells, transfer the cell culture from cell culture flasks to appropriate tubes. Centrifuge.
- Remove supernatant by decantation. Resuspend the cells in an appropriate volume of semi-complete medium (the volume is dependent on the system used for counting the cells).
- Count the cells.
- Resuspend the cell culture in semi-complete medium to a volume corresponding to a cell concentration of 2×10^5 cells /ml.
- Add 0.26 μl GM-CSF (150 $\mu\text{g}/\text{ml}$) per 1 ml of cell suspension to the cell culture.
- Incubate in cell culture flasks.

Cryopreservations

- To freeze the cells, transfer the cell culture from cell culture flasks to 50 ml tubes and centrifuge.
- Remove supernatant by decantation. Resuspend the cells in an appropriate volume of semi-complete medium.
- Count the cells. Resuspend the cell culture in semi-complete medium to a volume corresponding to a cell concentration of 14 million cells /ml.
- Add 0.26 μl GM-CSF (150 $\mu\text{g}/\text{ml}$) per 1 ml of cell suspension to the cell culture (final GM-CSF is 40 $\text{ng}/\mu\text{l}$)
- Prepare a solution of complete medium supplemented with 20% v/v DMSO.

- Transfer 0.5 ml of cell suspension to cryogenic vials (marked with cell bank identity, i.e. name of cells, date of freezing, type of cell bank stock).
- Add 0.5 ml of DMSO-supplemented complete medium to each of the cell-containing cryogenic vials, close the lids and invert the vials to mix.
- Immediately freeze the cells slowly in a temperature-controlled manner (-1°C/min down to -70 – -90°C).
- Vials are submerged into liquid nitrogen for long-term storage.

Preparing Flow cytometry samples

All washing steps are performed in Wash buffer. All centrifugations are performed at 300-315xg, 5 min, 2-8°C. All incubations are performed in dark at 2-8°C. Each lot of mAbs needs to be titrated to determine antibody concentration giving saturation (see Annex 2. Antibody titration).

Note:

- Removal of supernatant during preparation of flow cytometry samples is done by aspiration, e.g. by pipetting or by using a vacuum system, **not** by decantation.

Phenotypic Quality Control

The same day as performing a chemical stimulation, the cells are quality controlled with a phenotypic analysis. This is done to ensure cells are maintained in an inactivated state and to detect phenotypic drift.

Count cells and prepare 6 flow cytometry samples, with 2×10^5 cells in each sample.

- Wash the cells by adding ~1 ml Wash buffer and centrifuge at 300-315xg, 5 min, 2-8°C.
- Remove the supernatant by aspiration and repeat the washing step. Resuspend in 50 µl Wash buffer.
- Stain cells as indicated in Table 5 by adding titrated mAbs or viability stain to each sample.

Table 5. Antibodies and viability stain used in the Phenotypic Quality Control.

Sample 1	Isotype FITC	Isotype PE
Sample 2	CD86-FITC	CD54-PE
Sample 3	HLA-DR-FITC	CD80-PE
Sample 4	CD34-FITC	CD14-PE
Sample 5	CD1a-FITC	
Sample 6	Propidium Iodide (PI)	

- Incubate in dark at 2-8°C for ~15 min.
- Wash the cells by adding ~1 ml Wash buffer and centrifuge the samples at 300-315xg, 5 min, 2-8°C.
- Remove the supernatant by aspiration, resuspend in appropriate volume Wash buffer (using the mAbs and flow cytometer suggested in this protocol, 200 µl is used).

Analyze the samples on a flow cytometer according to manufacturer's instructions (for FACSVerse, use flow rate 60-120 µl/min) Record approximately and at least 10,000 events and analyze using the gating instructions below.

Analysis of Cell population

Exclude dead cells and cell debris by setting the “Cells” gate in the FSC-A/SSC-A scatter plot using Sample 1 (Isotype control), see Figure 2. Apply the “Cells” gate on Sample 2-6.




Figure 2. Instructions for setting the gate for the SenzaCell cell population.

Analysis of Phenotypic Quality Control markers

Show the gated “Cells” population in a PE/FITC scatter plot. Set quadrants for PE and FITC positive and negative cells using Sample 1 (isotype controls) as Figure 3A. Apply the quadrant from the isotype control sample in a PE/FITC scatter plot showing the “Cells” population of Sample 2-5 (mAb stained). Calculate and record the fraction of PE and FITC positive cells for each phenotypic marker, see example of Sample 4 below in Figure 3B, and compare with the accepted range in Table 6.




Figure 3.. Instructions for setting the quadrants for PE and FITC positive cells (A). Apply the preset gate and quadrants to record the fraction of positive cells for each phenotypic marker (B).

Analysis of Absolute viability (PI negative cells)

Set the gate for “Absolute viability”, in the PE/FITC scatter plot showing “All events” on Sample 1 (Figure 4A). Apply the “Absolute viability” gate on Sample 6 (PI stained), as in Figure 4B. Record the fraction of “Absolute viability” in % (PI negative cells) from Sample 6 and compare with the accepted range in Table 9.




Figure 4 Instructions for setting the gate for "Absolute viability" on the Isotype stained sample (A). Apply the gate "Absolute viability" on Sample 6 (B).

Phenotypic Quality Control acceptance criteria

The accepted range of phenotypic biomarker expression (Table 7) is based on observations made in the developing laboratory during assay development. Variations within these ranges are to be considered normal. However, if any of the biomarkers are out of the specified ranges, it is recommended that the cell batch is not used for cell stimulations at that timepoint. A representative example of a typical SenzaCell phenotype is presented in Annex 3. SenzaCell Phenotype. Note that the cell line is known to be heterogenous, and variations within the accepted ranges from given examples are expected.

Table 7. Accepted range for phenotypic markers of the SenzaCell cells.

Phenotypic biomarker	Accepted range of positive cells (%) *
CD86	10-40
CD54	>0
HLA-DR	>0
CD80	<10
CD34	>0
CD14	>0
CD1a	>0
Phenotypic biomarker	Accepted range of PI negative cells (%)
Absolute viability (PI negative cells)	≥84.5

Cell seeding for Test chemical stimulations

Cells are seeded for stimulation directly following a cell split, i.e. test chemical stimulations are to be scheduled to coincide with routine cell culture maintenance. This has been noted by the GARDskin assay developers to be an important factor. The cell stimulations are initiated when a stable cell culture is established i.e. when at least a duplication of the cells between cell passages is seen, and depending on the purpose of the cell stimulation, at specific cell passage ranges:

- For *cytotoxicity assessment*, cells at passage number **P4-P16** are used.
- For *GARD Main Stimulation*, cells at passage number **P6-P12** are used.

Test Chemical Stimulation Procedures

Handling of the test chemical

A chemical that is to be tested for sensitization in the GARDskin assay is referred to as a “test chemical”. The test chemical should be stored according to instructions from the supplier, in order to ensure its stability. Weighing of the test chemical can be performed prior to the day of cell stimulation if stored correctly and the stability of the substance can be ensured. Dissolved test chemical should be prepared fresh on the day of cell stimulation. Test chemicals should be dissolved in a compatible vehicle as appropriate stocks of target in-well concentration, in this document referred to as a Stock A concentration, depending on maximum possible solvent in-well concentration.

To prepare a *solid* test chemical, calculate the weight (see Note below about minimum weight of the scale) needed for an appropriate volume according to Equation 1. The test chemical is weighed into a pre-tared micro tube or appropriate vessel. Based on the weighed amount, calculate the exact volume of solvent needed to reach the c_T , according to Equation 1.

$$v = \frac{m \cdot p \cdot 0.01}{M \cdot c_T} \quad (\text{Equation 1})$$

Where

- V is the volume to be added in L
- m is the exact weight added to the tube in g
- M is the molecular weight of the test chemical in g/mol
- p is the purity of the test chemical in %
- c_T is the desired target concentration in mol/L

To prepare a *liquid* test chemical, use Equation 2 to calculate a dilution factor and calculate the volume of the test chemical and solvent needed for an appropriate volume of the test chemical of Stock A. Dilute the stock by the dilution factor into a 1.5 ml micro tube in the appropriate solvent.

$$df = \frac{c_S}{c_T} \quad (\text{Equation 2})$$

Where

- df is the dilution factor
- c_S is the concentration of the stock in mol/L
- c_T is the desired target concentration in mol/L

Note:

- A test chemical is preferably to be defined by a known molecular weight, as appropriate GARD Input Concentration are defined by molar concentrations. However for complex mixtures of unknown composition, substances of unknown or variable composition, complex reaction products or biological materials this may not be possible. However, this limitation may be circumvented by a) the use of weight-based concentrations (ppm), since the vast majority of sensitizers are detected at <100 ppm of the GARDskin assay (Gradin et al.), or b) by an approximation of the mean molar weight of the complex mixture.
- If the molar concentration of a liquid test chemical is not given by the customer, calculate the molar concentration using the molecular weight, density and purity of the test chemical.
- If the density of a liquid test chemical is not available, weigh the test chemical.
- If the substance is too viscous for pipetting, weigh the test chemical.

Solubility assessment

Consider the physicochemical properties of the test chemical and use the GARDskin compatible solvents and their maximum in-well concentrations listed in Table 7.

DMSO or water is the preferred vehicle. If neither is suitable, e.g. the test chemical is not dissolved, unfavourable reactions are expected, or other reasons, then a different vehicle can be tested. Independent of vehicle, a 2-step dilution is tested first, if not suitable then the 1-step dilution is tested. In some cases, where the test chemical has high reactivity with water or other unfavourable reactions are expected, a 1-step dilution can be tested first. Solubility of Test Items should be evaluated at least by a visual inspection of the solution. The Test Item is considered dissolved when the solution is without precipitate or phase separations. Other solubility approaches may also be used to reach dissolved test chemicals.

2-step dilution

Stock A

- Prepare a Stock A with a c_T (desired target concentration of the test chemical) of 500 mM.
- Vortex extensively and apply heat (37 ± 5 °C), if necessary, to dissolve the test chemical.
- Sonication may also be used if necessary and if compatible with the substance.
 - If not soluble: The stock A 500 mM maximum limit is specified for a vehicle concentration of 0.1% in well, if a higher vehicle concentration can be used (see Table 4) then lower the concentration accordingly (see Table 2) to achieve an in well concentration of 500 μ M. If soluble, proceed to the last step below. If the Test Item is not soluble at a concentration corresponding to 500 μ M in-well in any vehicle, identify the highest possible concentration where the Test Item is soluble in a suitable vehicle. Use the highest soluble concentration as a starting point in the dilution range.
 - Perform a serial dilution in the chosen vehicle to get a suitable range of stock A concentrations. Vortex well between each step. An example dilution scheme when working with a substance freely soluble at 500 mM in vehicle and at 500 μ M in medium is seen in Figure 5 on page 21.

Table 8 . The GARDskin compatible vehicles and the highest in-well concentrations. In the validation study only water and DMSO at 0.1 % were used. Other vehicles and vehicle concentrations have been verified to work in the GARDskin assay by the developing laboratory.

Vehicle*	Max in-well concentration (%)
DMSO	0.5
Water	0.1
Acetone	0.1
DMF	0.1
EtOH	0.1
Glycerol	1
Isopropanol	0.25
DMF:Glycerol 4:1 (v/v %)	0.25
Complete medium	100

*All the vehicles should be of high grade (purity $\geq 99.5\%$).

Serial dilution

- Perform a serial dilution in the chosen solvent to get a range of Stock A concentrations (see example in Figure 5). Vortex well between each dilution step. Extensive vortex and heat ($37^{\circ}\text{C}\pm 5^{\circ}\text{C}$) can be applied.
- From Stock A, prepare a range of Stock B concentrations by adding appropriate volume of Stock A to semi-complete medium. Extensive vortex and heat ($37^{\circ}\text{C}\pm 5^{\circ}\text{C}$) can be applied. If the substance is poorly soluble in semi-complete medium (typically identified through observation of precipitation in Stock B), the highest soluble concentration in semi-complete medium is used as the highest in the dilution range.
- In addition, prepare a Stock B concentration of chosen solvent (neg ctrl) in semi-complete medium to achieve the corresponding in-well concentration of the solvent.

1-step dilution

If a test chemical has a maximum solubility in Stock B which is lower than expected with respect to the maximum solubility in Stock A (i.e. displays solubility issues in Stock B), and is **found to be non-toxic**, i.e. the relative viability is above 95.4 %, the operator can try to increase the in-well concentration by bypassing the Stock B step as follows:

1. Prepare Stock A at the maximum possible concentration and prepare a serial dilution.
2. Add appropriate volume of Stock A dilution directly to complete medium.
3. Note the highest possible test chemical concentration in complete medium that does not have solubility issues i.e. completely solved.
4. If the highest possible concentration is higher than that in Stock B, prepare cells as described in 'cell seeding' in the sections below.

Note:

- Make sure to secure the lids before heating and vortexing the test chemical.
- A test chemical with solubility issues should be used from the highest possible concentration, down to 1 μM in-well concentration of the dilution range.

Cytotoxicity Assessment

The GARDskin Input concentration can be established for the test chemical in a *cytotoxicity assessment*. For an efficient workflow, multiple test chemicals can be assessed in each stimulation experiment with shared controls, unstim ctrl and neg ctrl, which are included in each *cytotoxicity assessment*. The cytotoxicity assessment results are used as guidance to find the appropriate stimulation concentration to be used during the Main Stimulations.

To determine the GARDskin Input concentration for a test chemical, cell stimulations are performed in a range of in-well concentrations as determined appropriate based on solubility assessment. For details about weighing and calculation, see section *Handling of the test chemical*.

Day 1

Cell seeding

The below cell seeding protocol is an example of the cell seeding when using a two-step dilution of the test chemical. If other dilutions schemes are used, use another seeding practice, but make sure the final cell concentration in-well after addition of the test chemical is 2×10^5 cells /ml.

- Transfer the cell culture from cell culture flasks to appropriate size tubes (e.g. 50 ml tubes). Centrifuge at 300-315xg for 5 min at 2-8°C.
- Remove the supernatant by decantation. Resuspend the cells in an appropriate volume of semi-complete medium.
- Count the cells.
- Resuspend the cell culture in semi-complete medium to a volume corresponding to a cell concentration of 2.2×10^5 cells /ml (final cell concentration in wells after addition of test chemical will be 2×10^5 cells /ml). Add 0.26 μ l GM-CSF (150 μ g/ml) per 1 ml of cell suspension to the cell culture.
- Use 24-well plates and seed **1.8 ml** of cell suspension into the number of wells needed (see Figure 5).

Cell stimulation

For each cell stimulation, independent of number of test chemicals and number of plates, it is recommended to include two wells of unstimulated cells (unstim ctrl) and one well with the vehicle (negative control) (Figure 5). If the test chemicals are dissolved using different solvents, each solvent should be present as negative controls. If several different concentrations of one solvent are used, it is sufficient to use the highest concentration as negative control in this experiment.

- **Test chemical:** add 200 μ l of Stock B to the 1.8 ml cell suspension seeded in 24-well plates for the dilution range of each test chemical. Mix well by pipetting. Final cell concentration in wells is 2×10^5 cells /ml.
- Example for 1-step dilution: If the final vehicle concentration should be 0.1% use the following method:
 - o Add 2 μ l of each Stock A dilution that did not have solubility issues in complete medium to wells containing 1.8 ml cells.
 - o Add 198 μ l semi-complete medium to the same wells and mix by pipetting.
- **Negative control:** add 200 μ l of relevant solvent Stock B to the 1.8 ml cell suspension. Example for 1-step dilution: If the final concentration should be 0.1% solvent use the following method:
 - o Add 2 μ l of solvent to wells containing 1.8 ml cells.
 - o Add 198 μ l semi-complete medium to the same wells and mix by pipetting.
- **Unstimulated control:** add 200 μ l of semi-complete medium to the 1.8 ml cell suspension, to achieve an in-well cell concentration and total volume equivalent to test chemical treated samples.
- Incubate for 24 h \pm 0.5 h at 37°C \pm 1°C and 5% \pm 0.5% CO₂.

Example of preparation of a serial dilution and a cytotoxicity assessment stimulation for one test chemical is shown in Figure 5.




Figure 5. Example of a schematic description of preparation of dilution series of one test chemical for cytotoxicity assessment.

Day 2

Cell harvest and sample preparation

After 24 h \pm 0.5 h incubation, harvest and prepare:

- Duplicate flow cytometry samples of each test chemical stimulation for the dilution range.
- Duplicate flow cytometry samples of negative control.
- Four flow cytometry samples of unstimulated cells (unstim ctrl).
- Mix the cell cultures in each well by carefully pipetting up and down and split into duplicate samples, \sim 1000 μl to each replicate. For the unstim ctrl four flow cytometry samples should be generated; two for staining with PI and two for non-staining. It is recommended to clearly mark the two samples which should not be stained. If only one unstim ctrl well has been used, split this cell culture into four flow cytometry samples.
- Wash the cells by adding \sim 1 ml Wash buffer and centrifuge.
- Remove the supernatant by aspiration, resuspend in \sim 1 ml wash buffer and centrifuge.
- Prepare a staining solution (enough for 50 μl for each flow cytometry sample) of 50 μl Wash buffer and 1 μl Propidium Iodide (PI).
- Resuspend each sample in 50 μl of the staining solution. **Note:** Leave 2 (out of 4) samples with unstimulated cells unstained, resuspend them in 50 μl Wash buffer instead.
- Incubate in dark at 2-8 $^{\circ}\text{C}$ for \sim 10 min.
- Wash the cells by adding \sim 1 ml Wash buffer and centrifuge. Resuspend in a suitable amount of Wash buffer (using FACSVerse: resuspend in 200 μl).

- Analyze the samples on a flow cytometer according to manufacturer's instructions (FACSVerse: flow rate 60-120 $\mu\text{l}/\text{min}$). Record approximately and at least 10,000 events and analyze using the gating instructions below.

Analysis of Cell population

Exclude dead cells and cell debris by setting the "Cells" gate in the FSC-A/SSC-A scatter plot using the unstimulated unstained sample, see Figure 6. Apply the "Cells" gate on all PI stained samples and record fraction of "Cells".

Note:

- The "Cells" population is not used for further analysis but is used to keep track of the placement of the cell population in the FSC-A/SSC-A scatter plot, see Annex 4. Cell population for a common pitfall. A low "Cells" population can give a false percentage "Absolute viability".




Figure 6. Instructions for setting the gate for the SenzaCells.

Analysis of Absolute viability (PI negative)

Use the unstimulated unstained sample to set the gate for "Absolute viability", in the PE/FITC scatter plot showing "All events" (Figure 7A). Apply the preset "Absolute viability" gate on all PI stained samples as in Figure 7B. Record the fraction of "Absolute viability" in % (PI negative cells).




Figure 7 Instructions for setting the gate for "Absolute viability" (A). Apply the gate "Absolute viability" on all PI stained samples (B).

Once the fraction of Absolute viability in % for the entire dilution range of a test chemical has been recorded, the Relative viability for each sample is calculated according to Equation 3. For each concentration of the dilution range, calculate the mean value of the duplicate samples.

$$Rv = \frac{V_s}{V_c} \cdot 100 \quad (\text{Equation 3})$$

Where

Rv is the **Relative** viability of the sample in %

V_s is the **Absolute** viability of the sample in %

V_c is the mean **Absolute** viability of the two unstimulated PI stained control samples in %

Note:

- The controls should pass the viability acceptance criteria, see Table 4, unstim ctrl: **Absolute** viability ≥84.5% and neg ctrl: **Relative** viability ≥95.5%.

The **GARDskin input concentration** of a test chemical is decided as following:

1. A test chemical that induces cytotoxicity should be used for *GARDskin Main Stimulation* at the concentration that induces 90% Relative viability (Rv90), where an acceptance criterion for each sample is a Relative viability of 84.5%-95.4%. If multiple concentrations fulfill the acceptance criterion, the concentration that yields the Relative viability **closest to 90%** is chosen as the GARDskin input concentration.
2. If the Relative viability decreases from ≥95.5% to <84.5% between two data points within the dilution range, the cytotoxicity assessment is repeated with a number of data points within the relevant concentration range. Interpolation between data points is not recommended, as linearity cannot be assumed.
3. A test chemical that is not cytotoxic (Relative viability ≥95.5%) is used for *GARD Main Stimulation* at 500 μM or highest soluble concentration.

GARDskin Main Stimulation

The purpose of the Main Stimulation is to generate three biological replicates of cell cultures exposed to the test chemical. The three should have been stimulated with the same test chemical concentration. However, in Main Stimulations several concentrations of the same test chemical may be used. For example time restrictions may necessitate the cytotoxicity assessment to coincide with the Main Stimulations.

The Main Stimulations should be performed three times with individual preparations of the test chemicals and controls. Individual cell batches (see section Abbreviations and definitions) should be used, to achieve three biological replicate samples (see Table 4 for details of the controls). The chemical concentration to be used in this step could be determined during the cytotoxicity assessment or directly in the main stimulation step where several concentrations may be assayed. The three *GARDskin Main Stimulations* can either be run in parallel or sequentially. If several test chemicals are to be analyzed in the same GARDskin experiment, the controls can be shared. Prior to starting a Main Stimulation, it is recommended to assign all Test and Reference Items to be tested a Study-unique Sample ID according to standard procedure set up in the testing laboratory. Consult permissible characters for the endpoint measurement CDF-files (NanoString technologies) and avoid incompatible characters in the Sample ID. Use the Sample ID for the Main Stimulation and for further procedures.

In Figure 8, a schematic example of a stimulation experiment with 8 test chemicals and the 3 controls are visualized, including one extra cell culture well with unstimulated controls.

Day 1

Preparation of test chemical and controls

- Prepare the test chemical in the same way as determined during solubility assessment and cytotoxicity assessment.
- Take into account the increased volume, 4 ml final cell suspension, for the stimulation.
- In addition, prepare the positive and negative control to achieve appropriate in-well concentration.

Cell seeding

Below cell seeding protocol is an example of the cell seeding when using a two-step dilution of the test chemical. If another dilution scheme was determined during solubility and cytotoxicity assessment, other cell seeding practices apply. No matter how the test chemical is diluted, the final cell concentration in wells after addition of test chemical should be 2×10^5 cells /ml.

- Transfer the cell culture from cell culture flasks to appropriate tubes. Centrifuge at 300-315xg for 5 min at 2-8°C.
- Remove supernatant by decantation. Resuspend the cells in an appropriate volume of semi-complete medium.
- Count the cells.
- Resuspend the cell culture in semi-complete medium to a volume corresponding to a cell concentration of 2.2×10^5 cells /ml (final cell concentration in wells after addition of test chemical will be 2×10^5 cells /ml). Add 0.26 µl GM-CSF (150 µg/ml) per 1 ml of cell suspension to the cell culture.
- Use 12-well plates and seed 3.6 ml of cell suspension into the number of wells needed for test chemicals and controls (see Figure 8).

Cell stimulations

For each cell stimulation experiment, it is recommended to include two cell culture wells of unstimulated cells (unstim ctrl). Below is, for convenience, described a test chemical standard dilution. Since the test chemical preparation can differ this is only a recommendation and guidance.

- **Test chemical:** add 400 µl of Stock B to the 3.6 ml cell suspension seeded in 12-well plates. Mix well by carefully pipetting up and down. Final cell concentration in well is 2×10^5 cells /ml.
If the 1-step dilution method is used (example for 0.1 % vehicle in-well):
 - o Add 4 µl of each Stock A to wells containing 3.6 ml cell suspension.
 - o Add 396 µl semi-complete medium to the same wells and carefully pipette up and down to mix.
- **Positive control:** add 400 µl of Stock B to the 3.6 ml cell suspension.
- **Negative control:** add 400 µl of Stock B to the 3.6 ml cell suspension.
If the 1-step dilution method is used and the final solvent concentration should be 0.1 %:
 - o Add 4 µl of solvent to wells containing 3.6 ml cell suspension.
 - o Add 396 µl semi-complete medium to the same wells and carefully pipette up and down to mix.

- **Unstimulated control:** add 400 µl of semi-complete medium to the 3.6 ml cell suspension, to achieve an in-well cell concentration and total volume equivalent to test chemical treated samples.
- Incubate for 24 h±0.5 h at 37°C±1°C and 5%±0.5% CO₂.




Figure 8. Schematic visualisation of test chemical addition to a GARDskin Main Stimulation.

Day 2

Cell harvest for RNA sample preparation and flow cytometry

From each cell culture well, prepare samples for RNA preparation and flow cytometry analysis as described below.

Samples for RNA preparation:

- Resuspend the cells by carefully pipetting up and down, ~1 ml volume, swirl on the bottom.
- Harvest 2 – 3 ≥1 ml samples from each cell culture well into separate RNase free 1.5 - 2 ml tubes and keep the tubes at 2 – 8°C for now.
- Harvest all samples from all plates before continuing.

Samples for flow cytometry analysis:

- Harvest the remaining volume from each cell plate well and split into two flow cytometry samples (~500 µl each).
- For the unstim ctrl four samples should be generated. Two for PI staining and two for non-staining. If only one unstim ctrl well has been used, split this well into four flow cytometry samples.

Note: Recommendation: clearly mark the two samples which should not be stained.

- Store the flow cytometry samples at 2 – 8°C.
- Harvest all samples for flow cytometry from all plates before continuing.

Samples for RNA preparation:

- Centrifuge the samples for RNA preparation at 300-315×g for 5 min at 2-8°C.
- Remove supernatant carefully by aspiration. Optional: depending on format aspirate only up to 12 tubes at a time, to avoid long-term contact with air and degradation of the RNA.

- Quickly add 500 µl of TRIzol or equivalent reagent, as specified in the RNA isolation kit protocol, to each cell pellet.
- Homogenize cells by vortexing the samples for 10-20 sec.
- Homogenized samples can be stored short-term in RT for maximum 1 hour, long-term at -18°C or below (stable for one month), or -80°C ±10°C (stable for one year).

Note: Though several reagents can be used instead of TRIzol, for practical reasons, these samples are called TRIzol samples throughout this protocol.

- Samples for flow cytometry analysis: Wash the flow cytometry samples by adding ~1 ml Wash buffer and centrifuge at 300-315xg, 5 min, 2-8°C.
- Remove the supernatant by aspiration and repeat the washing step.
- Prepare a staining solution (enough for 50 µl for each flow cytometry sample) of 50 µl Wash buffer and 1 µl Propidium Iodide (PI).
- Resuspend each sample in 50 µl of the staining solution.

Note: the two unstimulated cells which are unstained are resuspended in 50 µl Wash buffer instead.

- Incubate in dark at 2-8°C for ~10 min.
- Wash the cells by adding approximately 1 ml Wash buffer and centrifuge at 300-315xg, 5 min, 2-8°C.
- Remove the supernatant by aspiration and resuspend in appropriate volume (FACSVerse: 200 µl) Wash buffer.
- Analyze the samples on a flow cytometer according to manufacturer's instructions (FACSVerse 60 – 120 µl/min). Record approximately and at least 10,000 events.

Analysis of Cell population

Analyze the flow cytometry results in the same way as for the cytotoxicity measurement and calculate the mean Relative viability for each substance using Equation 3.

The purpose of the propidium iodide (PI) stained samples in the Main Stimulation is a **Quality Control of the Relative viability** to ensure that the test chemical and controls passes the relative viability acceptance criteria.

Relative viability acceptance criteria

- The unstimulated control should have an Absolute viability of ≥84.5%.
- The negative control should have a Relative viability of ≥95.5%.
- The positive control should have a Relative viability of ≥84.5%-95.4%.
- Test chemicals that are expected to induce cytotoxicity based on the cytotoxicity assessment should have a Relative viability between 84.5%-95.4%.
- Test chemicals that are assayed at 500 µM, or at the highest soluble concentration, should have a Relative viability of ≥84.5%.

If a test chemical or control does not pass the Relative viability acceptance criteria, either the entire Main stimulation should be discarded or only the test substance replicate from a single Main Stimulation should be discarded according to the decision tree below.

Procedures at failed Relative viability Quality Control

- If the Quality Control criteria is not reached for the **unstimulated control** (Absolute viability ≥84.5%), the Main stimulation has failed and all samples are discarded. A new Main Stimulation should be performed.

- If the Quality Control criteria is not reached for **the positive control** (Relative viability 84.5% – 95.4%) the Main stimulation has failed and all samples are discarded. A new Main Stimulation should be performed.
- If the Quality Control criteria is not reached for **the negative control** (Relative viability $\geq 95.5\%$), the Main stimulation has failed and all samples are discarded. A new Main Stimulation should be performed.
- If the Quality Control criteria for the Relative viability is not reached for a **test chemical**, discard its generated samples from this stimulation. The controls for this Main stimulation may still be used for other test chemicals within the same Main stimulation, but are not used for this test chemical. Include the test chemical in a new Main Stimulation. If necessary, the Test chemical may be re-analysed for cytotoxicity as described above in *Cytotoxicity Assessment*, before performing a new Main Stimulation.

Part result:

For *each* test chemical and control, 3 replicate samples for RNA isolation with passed Relative viability Quality Control are required. The replicates are generated in three individual *GARDskin Main Stimulations* (*i.e. all samples for RNA isolation generated from test chemical stimulations and control stimulations have to pass the viability quality control criteria within the same Main stimulation*).

RNA isolation

The RNA isolation is recommended to be performed with maximum 24 samples at a time. One vial of the replicate TRIzol samples generated from each Main Stimulation well is usually sufficient to obtain enough RNA for the endpoint measurement. If, for any reason, the first isolated sample does not meet the acceptance criteria, a second sample is generated, either by isolating the final TRIzol sample, or, when 3 or more TRIzol samples were generated from one well, by pooling two TRIzol samples onto the same spin column. For the second RNA sample isolated, the sample ID is suffixed, distinguishing it from the first, failed, sample.

Total RNA is isolated from the TRIzol samples using a commercially available kit and reagents. Direct-zol RNA MiniPrep, Zymo Research, specified in Table 3, is used by the assay developers and is therefore recommended. In general, the manufacturer's instructions should be followed.

Prepare the buffers and follow the protocol in the instruction manual included in the recommended kit (Direct-zol RNA MiniPrep), but with the following adjustments:

- Thaw TRIzol samples on ice.
- DNase I treatment should not be performed.
- After the centrifugation with RNA Wash buffer, discard the flow-through of the RNA Wash Buffer (re-use the collection tube) and perform an *additional* 1 min centrifugation (10,000-16,000xg) to avoid RNA Wash Buffer residues in the eluate.
- Elute RNA by adding ≥ 25 μ l DNase/RNase free water directly to the column matrix and centrifuge at 10,000-16,000xg for 30 seconds (though not recommended in the DirectZol protocol, the RNA can be eluted with as little as 20 μ l of DNase/RNase free water). If necessary to increase RNA yield, it is possible to perform a double elution by loading the eluted RNA once again on the same column for a second centrifugation.
- The eluted RNA can be used immediately or stored at -70 – -90°C.

It is recommended that a small aliquot (2 µl) is stored separately or used immediately, for quantification and quality control purposes.

RNA quantification and quality control

Analyze the RNA from each sample using an Agilent Bioanalyzer, or an equivalent instrument. Follow protocols provided by the supplier. RNA concentration and quality should correspond to NanoString recommendations. During test method development and validation, a sample with an RNA Integrity Number (RIN) of 8.0 and above, as derived from the Agilent Bioanalyzer 2100, was considered a sample of high quality. Corresponding or otherwise equivalent RNA quality metrics may be used to assure high quality RNA.

Part result:

For each test chemical and control, three RNA samples, each identified with a unique Sample ID and passing the RNA Quality Control is generated.

GARDskin Endpoint measurement

The endpoint measurement of GARDskin is the mRNA quantification of the GARDskin prediction signature (GPS) using the NanoString nCounter system. In the provided instructions below the nCounter MAX system is used. If another equivalent nCounter platform is to be used, refer to the manufacturer's instructions for appropriate experimental procedure and discuss with method developer.

The custom made CodeSet (i.e. a set of oligonucleotide probes representing the genes of the prediction signature) has been developed by SenzaGen and NanoString. To place an order for a batch of the CodeSet, please contact SenzaGen AB. It is not possible to use different batches of CodeSets for analysis of the endpoint using the GDA application. The nCounter analysis is performed with maximum 12 RNA samples at a time (one cartridge and CodeSet).

Setting up a NanoString Hybridization assay

All hybridization reactions use a total RNA input of 100 ng. According to the protocol below, the sample is added to the reaction in a volume of 5 μ l. Thus, all samples are to be diluted to a concentration of 20 ng/ μ l.

General Probe Handling warning: During setup of the assay, do not vortex or pipette vigorously to mix as it may shear the Reporter CodeSet. Mixing should be done by flicking or inverting the tubes. Do not spin any faster than 1000 \times g for more than 30 seconds as this may spin the CodeSet out of solution. Be aware that the maximum speed of a minicentrifuge is usually $>1000 \times$ g and that the "pulse" option of a microcentrifuge quickly goes to $>1000 \times$ g.

- Heat a thermocycler to 65°C, with a lid temperature of 70°C. A time program can be used. **Note:** Program the thermocycler using 15 μ l volume, at temperatures stated above.
- Thaw Reporter CodeSet and Capture ProbeSet at room temperature and store on ice. Flick and spin down. **Note:** After thawing, inspect the tube of Reporter CodeSet to make sure no colored precipitate is present. If you see a colored precipitate, heat the entire tube to 75°C for ~10 min and cool at room temperature before using.
- Thaw RNA samples on ice. Flick and spin down.
- Dilute all RNA samples to a concentration of 20 ng/ μ l using RNase-free water. Label each tube with its Sample ID. Mix by flicking and inverting and spin down.
- Prepare a master mix by adding 70 μ l of the hybridization buffer (provided in the NanoString master kit) to the Reporter CodeSet. Carefully mix by flicking and inverting and spin down.
- Cut a 12-strip of hybridization-tubes (provided in the NanoString master kit) in half, if necessary to fit them into a centrifuge. Note that the hybridization-tube strip has an orientation from 1-12, shown by the indent after the 1st and 8th position. Mark the tubes with Sample ID.
- Distribute 8 μ l of the Reporter CodeSet master mix to each hybridization tube.
- Add 5 μ l of diluted RNA sample to each hybridization tube. Carefully mix by flicking and inverting and spin down.
- Add 2 μ l of the Capture ProbeSet to each hybridization tube. Close the tubes with plastic lids and carefully mix by flicking and inverting and spin down.
- Place the hybridization tubes in the preheated thermocycler and incubate for 24 \pm 0.5 h.

Setting up a nCounter Prep Station Run

- Remove the cartridge from storage at -20 \pm 4°C and equilibrate to RT for ~30 minutes before the seal is broken.

- Remove 2 reagent plates from storage at 4°C and equilibrate to RT for ~30 minutes. Centrifuge the reagent plates 2000xg for 2 min, remove plastic lids before installing in the Prep Station.
- Chose “Start processing” on the screen and select the high sensitivity mode.
- Install all components into the Prep Station according to instructions on-screen. All required plastic material is provided in the NanoString master kit “prep pack” (stored at RT).
- After 24h (± 0.5 h) hybridization, remove the hybridization tubes from the thermocycler and spin down. Lift the metal lid in the prepstation and place the tubes without lids in the position highlighted on-screen. Orient the tubes so that sample #1 is positioned to the left and position #12 is positioned to the right. Close the metal lid.
- Initiate the run by following on-screen instructions. It is recommended to supervise the prep station until it starts piercing the foil of the prep plates. If reagents or plastics are not inserted correctly the Prep Station will have detected it at this point and informed on the screen.
- Once the Prep Station protocol is finished, after approximately 3 h, carefully remove the cartridge and place it on a lab-tissue and seal it with a provided adhesive cover.
- Discard used material in the Prep Station.
- Proceed to *Analysis with Digital Analyzer*.

Analysis with Digital Analyzer

- Create a Cartridge Definition File (CDF) according to NanoString’s instructions, using the highest field of view count, and upload it to the Digital Analyzer. The file maps the reads from each lane in the cartridge to sample specific attributes. Upload the created CDF-file to the Digital Analyzer prior to starting the analysis .
- Upload the RLF file to the Digital analyzer if it has not been added previously prior to initiating the analysis.
- Place the cartridge in the Digital Analyzer chose “Add Cartridge” and follow instructions on screen.

Note: If the Digital Analyzer is quantifying another cartridge, press pause and follow the instructions to add the new cartridge.

- Follow the on-screen instructions to start the analysis, the analysis takes ~5 hours for a full chip and the instrument can run overnight.
- When the Digital Analyzer is finished, download your RCC files according to internal documentation (FTP, email, and USB transfer of data is available). If stored at 2-8°C, the cartridge can be rescanned at least once for at least two weeks.

RCC file Quality Control

The output from each sample analyzed in the nCounter will automatically be quality controlled with internal control probes included in the CodeSet.

For NanoString facilities not using the GARD Data Analysis Application (GDAA), each acquired RCC-file should be quality controlled to assure that the nCounter analysis has been successful. Samples that fail any of the below described critical quality criteria should not be used for further analysis in the GARD™ skin data analysis pipeline. The critical quality metrics are imaging quality, linearity of the spike-in RNA control probes, limit of detection (LOD), and binding density. The imaging quality is calculated as the ratio between the number Fields of Views (FOV) (predefined in the CDF file as 555) and the number of successfully counted FOVs. A ratio above 0.75 (>0.75) is required for a sample to pass the imaging quality control. The linearity of the positive spike-in controls is calculated using the positive control probes

(POS_A-E) and their known RNA concentrations. The acquired counts for the positive control probes and their respective concentrations should be logarithmized (log2) before calculating the R² value of a linear fit to the data points. An R² value above 0.95 (R²>0.95) is required for a sample to pass the linearity quality control. The LOD quality control uses all the negative controls (NEG_A-H) and the positive control E (POS_E). The LOD is defined as the mean counts of the negative control probes plus 2 standard deviations of the counts, see Equation 4.

$$LOD = \mu + 2 * \sigma \quad (\text{Equation 4})$$

Where μ and σ are the mean value and the standard deviation of the negative control probes' counts respectively. For a sample to pass the LOD quality control, the positive control probe POS_E must be above the estimated LOD (POS_E > LOD). The binding density is a measure of the number of probes observed per cartridge surface area during the gene expression acquisition in the Digital Analyzer. For a sample to pass the binding density quality control, the binding density must be above 0.05 and below 2.25 (0.05 < binding density < 2.25). For a summary of the critical quality control parameters, see Table 8.

Table 8. Summary of the critical RCC-file quality control parameters

Quality Metric	Critical parameter
Imaging Quality	>0.75
Linearity	>0.95
Limit of Detection	<POS_E
Binding Density	0.05 - 2.25

In addition to the above described critical quality control parameters, it is also recommended to count the number of endogenous probes with 0 observed gene counts. If any samples contain multiple endogenous genes with 0 observed counts, a plausible explanation could be that the cartridge was analyzed with the wrong version of the RLF file. If this is the case, the affected cartridge should be rescanned (within two weeks) with the correct version of the RLF file to maintain sample integrity, and the previously generated RCC files should be discarded. If no apparent cause could be identified for the presence of multiple endogenous genes with 0 observed gene counts, other troubleshooting actions should be taken (not described herein).

Part result:

For each Sample ID analyzed by NanoString, a NanoString raw data file (RCC-file) is generated.

Data Analysis with GDAA

GDAA is an application that facilitates the endpoint analysis and the RCC file Quality control. Contact the GDAA support at gdaa.support@senzagen.com to receive a user account and the appropriate web site to be able to log in to the latest version of the GDAA.

Each sample (except unstim ctrl) analyzed with the GARDskin prediction model by GDAA will be assigned a decision value (DV). The GARDskin classifications of test chemicals, negative and positive controls correspond to the mean decision values (DVs) of the three GARDskin biological replicates.

- For skin sensitizers the mean DV ≥ 0 (UN GHS category 1).
- For non-sensitizers the mean DV < 0 .

Important: The RCC files for the test chemical should be analyzed together with relevant control RCC files i.e. only control samples (unstim, neg and pos) that are stimulated in the same *GARD Main Stimulation* as the test chemical.

Log in to the GDAA in a web browser and upload RCC files and annotation files. Annotation files are prepared according to the instructions under "Information" in the GDAA. The GDAA will generate a report (PDF format) containing the classifications and quality control. To be able to verify the integrity of the report, a MD5 checksum is also provided and it is recommended to save this code.

Acceptance Criteria

For the final GARDskin prediction to be valid for a test chemical, the following acceptance criteria must be met by the test chemical:

Phenotypic Control

All cell stimulations should have been performed with a batch of SenzaCell cells that passed the acceptance criteria of the Phenotypic Quality Control (Table 9).

Table 9. Acceptance range for phenotypic markers of the SenzaCells.

Phenotypic biomarker	Acceptance range of positive cells (%)
CD86	10-40
CD54	>0
HLA-DR	>0
CD80	<10
CD34	>0
CD14	>0
CD1a	>0

GARD Main Stimulations

- a. Should have been performed in three *GARD Main Stimulations* where each *GARD Main Stimulation* generated one biological replicate of the test chemical.
- b. The generated replicates should have passed Relative viability Quality Control.

RNA and RCC file Quality Control

- a. At least three replicates should have passed the RNA Quality Control.
- b. At least three replicates should have passed the RCC file Quality Control.

GARDskin classifications

- a. The classification should be based on three test chemical replicates, each of which should be derived from a *GARD Main Stimulation* that **also** generated:
 - i. An unstimulated control that passed Absolute viability-, RNA- and RCC Quality Control.
 - ii. A positive control that passed Relative viability-, RNA- and RCC Quality Control.
 - iii. A negative control that passed Relative viability-, RNA- and RCC Quality Control.
- b. The positive control and negative control should be accurately classified as a skin sensitizer and a non-skin sensitizer, respectively, by the GARDskin prediction model, i.e. the mean decision value of the positive control should be ≥ 0 and the mean DV of the negative control should be < 0 .

If these criteria cannot be met for a test chemical, the necessary parts of the GARDskin assay needs to be repeated for the specific test chemical to meet the acceptance criteria.

Recording of data

Recommended documentation regarding test chemical information and essential results are listed below. For a full overview of what to include in a test report see TG xxx.

Test chemical

- Substance ID, i.e. name of the test chemical or a simplified code of your choice (avoid special characters) used throughout the study to identify the substance. For the controls use “pos ctrl”, “neg ctrl” and “unstim ctrl”
- Molecular weight, density and purity used in the calculations
- Selected solvent and solvent in-well concentration
- Max screened in-well concentration
- GARD Input Concentration (μM)
- Relative viability at 90% (Yes/No)

Test method conditions

- Cell passage number at cell stimulation
- Phenotypic Quality Control passed (Yes/No)

Cytotoxicity assessment

- Determined GARD Input concentration (μM).
- Relative (or Absolute for unstim ctrl) viability at GARD Input concentration.
- Relevant comments, e.g. if Stock B was bypassed with a 1-step dilution.

GARD Main Stimulations

- In-well test chemical concentration (μM).
- Relative (or Absolute for unstimulated control) viability (%).

RNA Isolation and RNA QC

For each RNA Sample:

- RNA concentration ($\text{ng}/\mu\text{l}$).
- RNA Integrity Number (RIN) or other relevant RNA quality unit.

Endpoint measurement

For each RNA sample analyzed with NanoString:

- NanoString Quality Control passed (Yes/No).

Data analysis

- Save the rendered report files of the predictions.

Results

- Number of replicates meeting acceptance criteria with respect to Relative viability, RNA QC and NanoString QC.
- GARDskin prediction results.

Test Report

See OECD TG 442E for full list of information to report.

Annexes

Annex 1. FBS assessment

Annex 2. Antibody titration

Annex 3. SenzaCell Phenotype

Annex 4. Cell population

Annex 1. FBS assessment

The SenzaCell cell line provided by SenzaGen AB to a licensed CRO are accompanied with a previously assessed lot of Fetal Bovine Serum (FBS). New lots of FBS must be assessed to assure that it can support the phenotypic characteristics of SenzaCells. This is done by culturing the cells for 10 passages (or longer) with a previously assessed lot of FBS in parallel with the new lot of FBS. Cell growth and phenotypic characteristics are analysed at regular intervals throughout the assessment. Thawing, maintenance of cells at every passage and Phenotypic QC are performed as previously described in this protocol.

The FBS assessment is performed as follows:

- Prepare semi-complete medium with both control FBS (Control semi-complete medium) and new FBS (New semi-complete medium).
- Thaw a vial of cells and seed in Control complete medium (Passage 0).
- The next day, count and seed half the cells in Control complete medium and the other half in New complete medium (Passage 1).
- For passages 2-3, maintain the cells and seed in 20 ml Control and New complete medium, respectively.
- The FBS assessment is initiated at passage 4. For passages 4-10, at each passage, count the cells, perform a Phenotypic QC and seed in 20 ml Control and New complete medium.
- After passage 10, if seven Phenotypic QC records have been obtained, the cells can be discarded.
- A prolonged evaluation may be necessary depending on the results.

Acceptance criteria

The new FBS lot must fulfil the following criteria for it to be accepted for use in a Study:

- Cells grown in the new FBS lot must pass all Phenotypic QC criteria in at least 5 of the 7 recordings.
- Cells grown in the new FBS lot must display a cell growth where the cell density is $\pm 30\%$ of the cell density in the control FBS in at least 5 of the 7 cell density recordings in passages 4-10.

Annex 2. Antibody titration

The antibodies to be used for phenotypic control of Senza Cells shall be titrated (Ab list in Annex 2 Table 1). The antibody titrations should be performed on SenzaCells everytime a different brand, product number or of if a new batch within the same brand and product number is used. The staining with antibodies should reach saturation while using as small amount of antibodies as possible. Cell preparation, staining incubation, and flow cytometry analyses are performed as in section Phenotypic Quality Control. Include the Isotype controls as gating reference during the titrations.

Annex 2 Table 1. Antibodies used in the Phenotypic Quality Control.

Sample	FITC conjugated antibodies	PE conjugated antibodies
Sample 1	Mouse polyclonal anti-IgG1-FITC	Mouse polyclonal anti-IgG1-PE
Sample 2	Mouse anti-human CD86-FITC	Mouse anti-human CD54-PE
Sample 3	Mouse anti-human HLA-DR-FITC	Mouse anti-human CD80-PE
Sample 4	Mouse anti-human CD34-FITC	Mouse anti-human CD14-PE
Sample 5	Mouse anti-human CD1a-FITC	None

The phenotypic biomarker antibody titration is performed as follows:

- Wash and prepare the cells for staining as in section Phenotypic Quality Control.
- Stain the cells with the antibodies by adding different volumes, for an example, see Annex 2 Table 2, including the concentration suggested by the supplier.
- If the volume for the antibody to be used for costaining is known (Annex 2 Table 1), costain the cells keeping the concentration of the costaining antibody constant. If both antibody concentrations are unknown, titrate both separately and then test the chosen concentrations costained, to detect any quenching effects.
- Choose the concentration of the new antibody at the saturation point.

The Isotype antibody titration is performed as follows:

When assessing one of the Isotype antibody controls, a titration range is performed to assess that there is no unspecific staining. If the isotype control antibody has been used previously (i.e. from same brand and same product number) it only has to be tested at the concentration previously used (e.g 2 µl) and one concentration above (see Annex 2 Table 2). Laboratory procedure is the same as for other antibody titrations described above.

Annex 2 Table 2. Example of titration volumes for Ab titration.


<i>Volume (μl) of old Ab and/or recommended volume according to Ab producer.</i>	0.5	1	1.5	2	2.5	3
Recommended titration volumes (μl) of new Ab.	0.25	0.5	0.5	1	1.5	2
	0.5	0.75	1	1.5	2	2.5
	0.75	1	1.5	2	2.5	3
	1	1.5	2	2.5	3	3.5
	1.25	2	2.5	3	3.5	4
	1.5*	2.5	3	3.5	4	4.5
	-	3*	3.5*	4*	4.5*	5*
Total volume of new Ab needed for titration (μl)	5.25	11.25	14	17.5	21	24.5

*If no saturation is seen, increase the Ab concentration.

The lowest concentration generating maximum fluorescence is selected. For representative figures of each antibody staining of SenzaCells see Annex 3. The GARDskin phenotypic control acceptance criteria is shown in Table 9.

Annex 3. SenzaCell Phenotype


An example of the SenzaCell phenotype is visualized in Annex 3 Figure 1.



Annex 3 Figure 1. PE/FITC scatter plots of mAb stained SenzaCell cells for phenotypic control.

Annex 4. Cell population

During analysis of the viable cell population it is important to keep track of the “Cells” population in the FSC(-A)/SSC(-A) scatter plot and the diagonal displacement of the “Absolute viability” population in the PE/FITC scatter plot. Annex 4 Figure 1 is visualizing the FSC(-A)/SSC(-A) scatter plot and PE/FITC scatter plot for cells stimulated with 400 μ M (A) and 100 μ M (B), respectively. A low fraction of “Cells” can give a false percentage of Absolute viability, see Annex 4 Figure 1A.



Annex 4 Figure 1. FSC/SSC and PE/FITC scatter plots of cells stimulated with a test chemical at 400 μ M (A) and at 100 μ M (B).

Bibliography

Scientific Publications

Johansson H., Gradin R., Johansson A., Adriaens E., Edwards A., Zuckerstätter V., Jerre A., Burleson F., Gehrke H., Roggen E. Validation of the GARD™skin assay for assessment of chemical skin sensitizers – ring trial results of predictive performance and reproducibility. *Toxicological Sciences*. May 17, 2019.

Gradin R., Lindstedt M., Johansson H. Batch Adjustment by Reference Alignment (BARA): Improved prediction performance in test sets with batch effects. *PLOSOne*, 2019.

Gradin, R., Forreryd, A., Mattson, U., Jerre, A. and Johansson, H. (2021). Quantitative Assessment of Sensitizing Potency using a dose-response adaptation of GARDskin. *Sci. Rep.* Manuscript in press.

Zeller K. S., Forreryd A., Lindberg T., Gradin R., Chawade A., and Lindstedt M. The GARD platform for potency assessment of skin sensitizing chemicals. *ALTEX Online* first published April 12, 2017.

Johansson H., Gradin R., Forreryd A., Agemark M., Zeller K., Johansson A., Larne O., van Vliet E., Borrebaeck C., Lindstedt M. Evaluation of the GARD assay in a blind Cosmetics Europe study. *ALTEX Online* first, 2017.

Forreryd A., Zeller K., Lindberg T., Johansson H., Lindstedt M. From genome-wide arrays to tailor-made biomarker readout – Progress towards routine analysis of skin sensitizing chemicals with GARD. *Toxicology In Vitro*, 2016.

Forreryd A., Johansson H., Albrekt A.S., Lindstedt M. Evaluation of high throughput gene expression platforms using GARD – an assay for prediction of skin sensitization based on a genomic biomarker signature. *BMC Genomics*, 2014.

Johansson H., Rydnert F., Kuehnl J., Schepky A., Borrebaeck C.A.K., Lindstedt M. GARD in-house validation – A proof of concept. *Tox Sci*, 2014.

Johansson H., Lindstedt M. Prediction of skin sensitizers using alternative methods to animal experimentation. *Basic Clin Pharmacol*, 2014.

Johansson H. GARD - Genomic Allergen Rapid Detection. A testing strategy for assessment of chemical sensitizers. *Academic Doctoral Thesis*, 2013.

Johansson H., Albrekt A.S., Borrebaeck C.A.K., Lindstedt M. The GARD assay for assessment of chemical skin sensitizers. *Toxicol In Vitro*, 2013.

Lindstedt M., Borrebaeck C.A.K. Pattern rules: biomarker signatures for sensitization as an alternative to animal testing. *Biomark Med*, 2011.

Johansson H., Albrekt A.S., Lindstedt M., Borrebaeck C.A.K. A genomic biomarker signature can predict skin sensitizers using a cell-based in vitro alternative to animal tests. *BMC Genomics*, 2011.