

final version 5 Novemb. 97	Standard Operation Procedure EpiDerm™ Phototoxicity Assay (model: Epi-200)	page 1 of 28
---	---	-------------------------------

ANNEX 3

Final SOP

final version 5 Novemb. 97	Standard Operation Procedure EpiDerm™ Phototoxicity Assay (model: Epi-200)	page 2 of 28
-------------------------------	---	-----------------

The following SOP is the final version used within the scope of the ECVAM project

"Evaluation of the Prevalidation Process"

subproject:

"Prevalidation of the EpiDerm™ Skin Phototoxicity Test"

- Phase III (blind trial) -

The first draft SOP based upon the Skin² test developed by ZEBET and Advanced Tissue Sciences for the full skin model in the EU/COLIPA joint project "In vitro Photoirritation". During phase I of the ECVAM prevalidation study ZEBET adopted the method to the EpiDerm™ technology and established a data base of 39 tests performed on 12 chemicals. Optimisation experiments revealed overnight exposure (18-24 hrs) with the test chemical and use of a light dose of 6 J/cm² the best design. Tests with chlorpromazine applied in H₂O, as well as in oil, and in a H₂O/oil emulsion showed that the assay is capable to handle formulations. A refined SOP (version 2, 30 May 1997) was distributed for comments to the labs participating in the ECVAM prevalidation study.


In a meeting in Berlin (4 September 1997) P&G, Beiersdorf and ZEBET agreed on several modifications of the SOP (version 2). The major modification was a new technique of topical application using paper pads. Since there was no experience with the pad technique additional identical experiments were performed in all of the 3 labs and in a phone conference held on 1 October 1997, P&G, BDF and ZEBET agreed on the final SOP to be used in phase III. Apart from minor changes in details and wording, compared to the 1st draft SOP, the main amendments of the final SOP comprise

- optionally, UV irradiation can be performed in 24 well plates on 0.3 mL medium instead of 6-well plates on 0.9 mL medium.
This change was made, since 0.3 mL medium is sufficient for supply of the tissues during 60 minutes irradiation. Thus, the weekly testing throughput can be increased.
- wherever possible, chemicals shall be applied as **solutions**, either in **oil** or in **H₂O**. If chemicals cannot be dissolved either in H₂O or in oil, they shall be applied as **suspensions in oil**.
- chemicals dissolved in **H₂O** are applied at **50 µL without using a pad**.
- chemicals dissolved (or suspended) in **oil** are applied at **20 µL using a pad** (Finn chamber disk, 8 mm Ø).
- reading of optical densities of formazan extracts is done with 570 nm (or equivalently 540 nm) **without** using a reference filter.
- a simplified Methods Documentation Sheet (MDS) is used
- a modified MS Excel data spreadsheet (P-SPREAD.XLS) is used

The following three laboratories have approved the final SOP and will perform testing in phase III according to the SOP:

- Dr Manfred Liebsch, ZEBET at the BgVV, Diedersdorfer Weg 1, D-12277 Berlin Germany
☎ +49-30-8412-2275 Fax +49-30-8412-2958 e-mail zebet@bgvv.de
- Dr. Frank Gerberick / Lynn Cruse, Procter & Gamble, 11810 East River Road ROSS, Ohio 45061, USA
☎ +1-513-627-2909 Fax +1-513-627-0400 e-mail gerberick.gf@pg.com
☎ +1-513-627-2909 Fax +1-513-627-0400 e-mail cruse.lw@pg.com

final version	Standard Operation Procedure	page 3
5 Novemb. 97	EpiDerm™ Phototoxicity Assay (model: Epi-200)	of 28

- Uwe Pfannenbecker, Beiersdorf AG, KSt 4232, Unnastraße 48, D-20245 Hamburg, Germany
 +49-40-4909-3916 **Fax** +49-40-4909-3589 **e-mail** 113157.26@compuserve.com

final version 5 Novemb. 97	Standard Operation Procedure EpiDerm™ Phototoxicity Assay (model: Epi-200)	page 4 of 28
---	---	-------------------------------

CONTENTS:

1 RATIONALE	6
2 NEED FOR THE ASSAY	6
3 BASIC PROCEDURE	6
4 MATERIALS	7
4.1 MATERIALS, <u>NOT PROVIDED</u> WITH THE KITS:	7
4.1.1 <i>Laboratory aids</i>	7
4.1.2 <i>Technical Equipment</i>	7
4.1.3 <i>UVA-vis Irradiation equipment</i>	8
4.1.4 <i>Solutions, Reagents</i>	8
4.1.5 <i>Computer software</i>	8
4.2 EPI-200 KIT COMPONENTS.....	8
4.3 MTT-100 ASSAY KIT COMPONENTS	9
5 METHODS.....	9
5.1 EXPIRATION AND KIT STORAGE	9
5.2 CALIBRATION OF THE SOLAR SIMULATOR	9
5.3 QUALITY CONTROLS.....	10
5.3.1 <i>UVA Sensitivity of the Epi-200 Tissues</i>	10
5.3.2 <i>Negative Control</i>	10
5.3.3 <i>Positive Control: Chlorpromazine (CPZ)</i>	10
5.3.4 <i>Maximum inter tissue viability difference of tissue couples</i>	11
5.4 TEST SAMPLE PREPARATION AND TEST CONCENTRATIONS.....	11
5.4.1 <i>Concentration series</i>	11
5.4.2 <i>Application of test sample</i>	12
5.5 EXPERIMENTAL PROCEDURE	12
5.6 DOCUMENTATION.....	15
5.6.1 <i>Method Documentation Sheet, MDS</i>	15
5.6.2 <i>Data Spreadsheet</i>	16
6 PREDICTION MODEL.....	16
7 REFERENCES	17
8 ANNEX A: PROCEDURE STEPS.....	18
9 ANNEX B: METHODS DOCUMENTATION SHEET	19
10 ANNEX C: EXCEL SPREADSHEET	21
11 ANNEX D: POSITIVE REFERENCE DATA.....	23
12 ANNEX E: NEGATIVE REFERENCE DATA	24
13 ANNEX F: EPI-200 UVA-SENSITIVITY	25
14 ANNEX G IRRADIANCE SPECTRUM OF THE SUN SIMULATOR.....	27

final version 5 Novemb. 97	Standard Operation Procedure	page 5 of 28
	EpiDerm™ Phototoxicity Assay (model: Epi-200)	

final version	Standard Operation Procedure	page 6
5 Novemb. 97	EpiDerm™ Phototoxicity Assay (model: Epi-200)	of 28

1 RATIONALE

Phototoxicity (photoirritation) is here defined as acute toxic response that is elicited after the first exposure of skin to certain chemicals and subsequent exposure to light, or that is induced similarly by skin irradiation after systemic administration of a chemical substance.

The present assay is designed to detect the phototoxic potential of a chemical by using a three dimensional human epidermis model*. Since the assay allows application of test materials to the air exposed surface (stratum corneum), it mimics the *in vivo* situation and thus may allow to predict phototoxic potency of test materials applied in usage concentrations. The test is based upon a comparison of the cytotoxicity of a chemical when tested with and without additional exposure to a non toxic dose of UVA+visible light. Cytotoxicity is expressed as reduction of mitochondrial conversion of MTT to formazan¹, determined one day after chemical treatment and UVA exposure.

* MatTek's EpiDerm System^{2, 3} consists of normal, human-derived epidermal keratinocytes which have been cultured to form a multilayered, highly differentiated model of the human epidermis. It consists of organised basal, spinous and granular layers, and a multi-layered stratum corneum containing intercellular lamellar lipid layers arranged in patterns analogous to those found *in vivo*. The EpiDerm tissues (surface 0.6 cm²) are cultured on specially prepared cell culture inserts (Millicells®, 10 mm Ø) and shipped worldwide as kits, containing 24 tissues on shipping agarose.

2 NEED FOR THE ASSAY

It has been shown in a joint EU/COLIPA validation project^{6, 7}, that the phototoxic potential of chemicals can be correctly predicted by using cell culture monolayers in a specially designed cytotoxicity assay, the 3T3-NRU-phototoxicity test. Since the phototoxic potential of a chemical predicted using a cellular system may not be relevant when topically applied to the skin at low concentrations (e.g. in a formulation) there is a need for adjunct tests, which allow for the assessment of safe usage concentrations on a dose per area basis before testing them in humans. Reconstituted skin models and epidermis models have shown to be able to predict both, photoirritancy^{4, 5, 8}, as well as the photoprotective action of sunscreens⁵. In addition, skin models can handle formulations (e.g. emulsions, suspensions) which the 3T3 test cannot handle. Thus, in a testing strategy which is based purely on *in vitro* tests, there is a need to combine the basic 3T3 NRU PT with other *in vitro* tests, which may allow to assess safety or phototoxic potency of formulations. In addition, a phototoxicity test involving a human skin model may be useful for risk benefit analysis of dermal pharmaceuticals.

3 BASIC PROCEDURE

On day of receipt (e.g. Tuesday afternoon) EpiDerm™ tissues are stored over night in a refrigerator. Next day, at least one hour before starting the assay, tissues are transferred to 6-well plates with assay medium and the medium is exchanged. Then, 5 concentrations of the test material (dissolved in H₂O or oil or suspended in oil) are topically applied onto 2 tissues per concentration (i.e. 1 vehicle control + 5 concentrations = 12 tissues). A second set of 12 tissues is treated identically. Plates are incubated over night. Next day, one set of tissues is exposed to 6 J/cm² UVA (+UVA **part of the test**) and the other set is kept in the dark for the same period (-UVA **part of the test**). Tissues are then rinsed with PBS to remove test material, transferred to new 6 well plates with fresh medium and incubated over night. Next day, assay medium is replaced by MTT-medium and tissues are incubated for 3 hours with MTT. Tissues are then rinsed with PBS, and the formazan is extracted with

final version 5 Novemb. 97	Standard Operation Procedure EpiDerm™ Phototoxicity Assay (model: Epi-200)	page 7 of 28
---	---	-------------------------------

Isopropanol. Optical density is determined at 540/570 nm in a plate spectrophotometer and cell viability is calculated for each tissue as % of the corresponding vehicle control either irradiated or unirradiated.

4 MATERIALS

4.1 *Materials, not provided with the Kits:*

4.1.1 *Laboratory aids*

Sterile, blunt-edged forceps	<i>for transferring inserts</i>
6-well tissue culture plates (in addition to those provided)	<i>If, instead of replacing media inserts are transferred to new plates with media</i>
24-well culture plates (in addition to those provided)	<i>for UVA irradiation, if not performed in 6-well plates</i>
96-well plates (flat-bottom)	<i>for OD reading in plate spectrophotometer</i>
Sterile disposable pipette tips	
5 Beakers á 50 ml	
sterile, capped glass or plastic test tubes	<i>for preparing the concentration series</i>
Finn chamber filter pads, Ø 8 mm, sterilised HERMAL, Scholtzstr. 3, D-21465 Reinbek, Purchase Order No.: D 9503	<i>for application of test materials dissolved or suspended in oil</i>
Repeat pipetter (2 ml)	<i>for adding the extractant solution</i>
Positive displacement pipettes (20 µL, 50µl)	<i>for application of viscous test materials</i>
adjustable Pipet (100 µl)	<i>for pipetting the concentration series</i>
adjustable Pipet (200 µl)	<i>for pipetting the concentration series</i>
adjustable Pipet (1000 uL)	<i>for medium change</i>

4.1.2 *Technical Equipment*

Bunsen burner or autoclave	<i>for sterilising forceps</i>
37 °C humidified incubator with 5% CO ₂ ,	<i>for incubating tissues prior to and during assays</i>
vacuum source/trap	<i>for aspirating solutions</i>
laminar flow hood	<i>for transferring tissues under sterile conditions and for application of test materials.</i>
37 °C water bath	<i>for warming up Assay Medium, PBS etc.</i>
Laboratory balance	<i>for preparing concentration series</i>
96-well Spectrophotometer (Plate-Reader) equipped with filter 570 nm or 540 nm	<i>for reading optical density at 570 nm,</i>
Shaker for cell culture plates	<i>for extraction of formazan</i>
Laboratory centrifuge 1500 x g	<i>for centrifugation of MTT medium</i>

final version 5 Novemb. 97	Standard Operation Procedure EpiDerm™ Phototoxicity Assay (model: Epi-200)	page 8 of 28
---	---	-------------------------------

Vortex mixer

for keeping test suspensions homogeneous during preparation of the concentration series

electric homogeniser

for preparing test chemical suspensions in oil

4.1.3 UVA-vis Irradiation equipment

UV-sun simulator
type SOL 500 or SOL 3

*Dr. K. Hönle GmbH, Frauenhoferstr. 5,
D-82152 Martinsried, Germany*

Contact: Dr. G. Schmid

☎: +49-89-856 08-0

Fax: +49-89-856 08-48

Any appropriate, adjustable and stable tripod

For the fixation of the SOL 500

UVA-meter, type No. 37, Dr. Hönle

For everyday check of calibration

UVA-meter, type No. 37, Dr. Hönle

Use only as a reference in case of unexpected readings with the everyday radiometer

Filter, type H1, Dr. Hönle

Use to cut-off emitted UVB

4.1.4 Solutions, Reagents

Sesame oil (USP/EP/DAB grade)
source: pharmacy

Vehicle for test materials

H₂O Aqua Pur (Millipore®), or distilled H₂O

Vehicle for test materials

PBS with Ca⁺⁺ and Mg⁺⁺: ~ 500ml per test
(e.g. Gibco # 14040)

For rinsing-off test materials after irradiation

4.1.5 Computer software

(MS Windows) software for Plate
spectrophotometer

software must be able to export data

MS Excel 5.0

For and data analysis in the Data Spreadsheet

4.2 Epi-200 Kit Components

Examine all kit components for integrity. If there is a concern call MatTek Corporation immediately (Mitch Klausner, ☎ +1-508-881-6771, Fax +1-508-879-1532).

1	Sealed 24-well plate	<i>Contains 24 inserts with tissues on agarose</i>
2	24-well plates	<i>Use for MTT assay and formazan extraction</i>
4	6-well plates	<i>Use for pre-incubation (and assay)</i>
1 bottle	Serum-Free Assay Medium	<i>DMEM-based medium</i>
1 bottle	Maintenance Medium	Do not use in the present assay
1 bottle	PBS Rinse Solution (100 mL)	<i>Use for rinsing the inserts in MTT assay</i>
1 vial	1% Triton X-100 Solution (10 mL)	<i>Skin irritant reference chemical</i> Do not use in present assay
1	MTT Assay Protocol	<i>steps are included in the present SOP</i>

final version 5 Novemb. 97	Standard Operation Procedure EpiDerm™ Phototoxicity Assay (model: Epi-200)	page 9 of 28
-------------------------------	---	-----------------

4.3 MTT-100 Assay Kit Components

1 vial, 2 ml	MTT concentrate	
1 vial, 8 ml	MTT diluent	For diluting MTT concentrate
1 bottle, 60 mL	Extractant Solution (Isopropanol)	For extraction of formazan crystals

5 METHODS

5.1 Expiration and Kit Storage

Epi-200 kits are shipped from Boston on Monday; make sure that they are arriving in the laboratory on Tuesday. Upon receipt of the EpiDerm tissues, place the sealed 24 well plates and the assay medium into the refrigerator (4°C). Place the MTT concentrate containing vial in the freezer (- 20°C) and the MTT diluent in the refrigerator (4°C).

part #	description	conditions	shelf life
EPI-200	EpiDerm™ cultures	refrigerate (4°C)	until Friday
EPI-100	assay medium	refrigerate (4°C)	7 days
MTT-099	MTT diluent	refrigerate (4°C)	7 days
MTT-100	MTT concentrate	freeze (- 20°C)	2 month

Record lot numbers of all components and transfer lot/production label on sealed tray onto the Methods Documentation Sheet (MDS see ANNEX B).

Note: Since testing starts on Wednesday, irradiation on Thursday, MTT assay on Friday, do not order more Epi-200 kits per week than can be dosed or irradiated on one day, respectively. This does not hold for US labs: if they receive kits on Tuesday before 12:00 the test can be started.

5.2 Calibration of the Solar Simulator

Note: New metal halide burners should be burned for ~100 hrs prior to first use to achieve a stable emittance. According to Dr Hoenle the burner has a shelflife (in which the spectrum is stable) of at least 800 hrs. Recording of lamp usage hours is, therefore, recommended. Extended use is only acceptable if the emitted energy spectrum can be checked.

1. Mount the SOL 500 / SOL 3 lamp, equipped with a H1-filter, on any appropriate stable tripod allowing fine-adjustment of the exposure distance.
2. Adjust SOL 500 / SOL 3 to a distance of about 60 cm.
3. Switch the Lamp on, wait at least 15 minutes and measure irradiance **through the lid of a cell culture plate** using the calibrated UV radiometer (type 37, Dr. Hönle), equipped with an UVA-sensor of the same serial number.
4. Adjust distance of SOL 500 / SOL 3 to achieve a UVA irradiance of **1.7 mW/cm²** (The resulting dose will be 1 J/cm² per 10 min. exposure time)
5. According to the number of plates to be exposed concurrently, check the exposure area for equal distribution of irradiance:

A range of **1.6 - 1.8 mW/cm²** is acceptable. **Important:** A maximum difference of **1.5 and 1.9 mW/cm²** can be accepted, if positions of the plates with low and high

final version 5 Novemb. 97	Standard Operation Procedure EpiDerm™ Phototoxicity Assay (model: Epi-200)	page 10 of 28
-------------------------------	---	------------------

irradiance are changed after half time of the irradiation (30 minutes) is reached (like chess castling).

Calibration of the SOL 500 / SOL 3 shall be checked as described above each time before performing a phototoxicity assay. In case measurements with the UV radiometer reveal unexpected results, either the metal halide burner may have reached the end of its shelflife, or the radiometer is de-calibrated due to various reasons. In this case, a second reference radiometer of the same type and calibration, which has not been handled every day and kept in the dark shall be used for cross check.

5.3 Quality Controls

5.3.1 UVA Sensitivity of the Epi-200 Tissues

Note: A UVA sensitivity experiment should be performed once the test is newly set up in a laboratory. If UVA sensitivity of the tissues is within the acceptance range this type of experiment should be repeated in greater intervals (e.g. once every 6 months).

1. Incubate 24 tissues (37°C, 5 % CO₂) for **at least 1 hr** in 6-well plates with 0.9 ml assay medium /well to allow release of metabolites and debris accumulated during the shipment.
2. Adjust irradiance of the SOL 500 / SOL 3 to **1.7 mW/cm²** (measure through plate lid !)
3. For UV irradiation, transfer 21 tissues to a 24-well plate filled with 0.3 mL assay medium per well.
4. Prepare a 24 well plate with 0.3 ml assay medium per well and transfer the 3 tissues serving as non-irradiated control. Place this plate in a dark box at room temperature.
5. Start irradiation of the 21 tissues **through the lid of the plate**. Use a fan to prevent H₂O condensation under the lid. Every 30 minutes (= 3 J/cm²) transfer 3 tissues from the irradiation site to the dark box. The resulting dose series is **3, 6, 9, 12, 15, 18, 21 J/cm²**.
6. Incubate tissues over night (18 - 24 hrs.) at 37°C, 5 % CO₂, 90 % humidity
7. Determine tissue viability according to **5.5**. Compared to the non irradiated tissues (100 % viability) up to 6 J/cm² (= 60 minutes) there shall be no reduction of viability exceeding 20%. The historical ID_{50 UVA} is in the range of ~12 - 18 J/cm² (see **ANNEX F**).

5.3.2 Negative Control

The **absolute OD** of the negative control tissues in the MTT-test (see **5.5**) is an indicator of tissue viability obtained in the testing laboratory after shipping procedure and under specific conditions of the assay.

*Tissue viability is meeting the acceptance criterion if the mean OD of the two negative control tissues (determined without reference filter) is **OD ≥ 0.8**.*

5.3.3 Positive Control: Chlorpromazine (CPZ)

For the present study, it is not necessary to include a positive control into each phototoxicity test as this reduces the number of concentrations of the test chemical. When the assay is newly established perform a full experiment with five concentrations of **Chlorpromazine** (dissolved in H₂O) ranging from 0.001% up to 0.1%. Repeat this test on a regular basis.

final version 5 Novemb. 97	Standard Operation Procedure EpiDerm™ Phototoxicity Assay (model: Epi-200)	page 11 of 28
-------------------------------	---	------------------

*A dose dependent reduction of cell viability occurring only in the UVA-irradiated tissues, shall be observed between 0.00316% and 0.0316% (see **ANNEX D**).*

Note: If, in other studies, CPZ shall be included in each assay, use 0.316%.

5.3.4 Maximum inter tissue viability difference of tissue couples

The new spreadsheet calculates differences in viability between tissue couples that are treated identically. According to the historical data base of ZEBET the mean difference between untreated tissue duplicates is 9% ± 7% (S.D.).

A difference > 30% (i.e. exceeding the 99% confidence interval) between two tissues treated identically should be regarded as a rejection criterion, and re-testing of the chemical is recommended if the resulting viability is near to the classification cut-off.

5.4 Test Sample Preparation and Test Concentrations

According to their solubility, chemicals are applied either as solution in **water**, or as solution or suspension in **sesame oil**.

DESCRIPTIVE TERM	RANGE OF SOLUBILITY	% (w/v)	CATEGORY
very soluble	> 1000 mg/ml	>100.00	1
freely soluble	> 100 mg/ml - 1000 mg/ml	>10.00	2
soluble	> 30 mg/ml - 100 mg/ml	>3.00	3
sparingly soluble	> 10 mg/ml - 30 mg/ml	>1.00	4
slightly soluble	> 1 mg/ml - 10 mg/ml	>0.10	5
very slightly soluble	> 0.1 mg/ml - 1 mg/ml	>0.01	6
practically insoluble	0.1 mg/ml and lower	<0.01	7

As a basic recommendation, poorly water soluble test materials (category 5-7) should be tested dissolved or suspended in sesame oil. Water soluble test materials (category 1-4) shall be tested dissolved (!) in water.

If suspensions are tested use appropriate techniques for preparing e.g. a homogeniser, Vortex or Sonicator).

5.4.1 Concentration series

Prepare five concentrations of the test material. Where possible, the highest concentration of a test material should show cytotoxicity in non-irradiated tissues. Since many test chemicals are likely to absorb UV they can act as UV-filter. Therefore, the highest test concentration should **not exceed 10%**.

If there is no information on skin toxicity of the test material, start with the following concentration series:

vehicle	% (w/v)	% (w/v)	% (w/v)	% (w/v)	% (w/v)
oil	10	3.16	1	0.316	0.1
water	1	0.316	0.1	0.0316	0.01

final version	Standard Operation Procedure	page 12
5 Novemb. 97	EpiDerm™ Phototoxicity Assay (model: Epi-200)	of 28

Note: According to ZEBET's experience, the series for oil fits for many test materials. Materials solved in water pass the stratum corneum more quickly. If they are, in addition, highly cytotoxic the concentration series may have to be shifted to a lower range in a second experiment.

5.4.2 Application of test sample

Note: To be as close as possible to the human *in vivo* test and to achieve a constant application area, a patch technique was developed using 8 mm Ø paper disks (normally used in Finn chambers for human patch tests). However, experiments performed at P&G, Beiersdorf AG and ZEBET revealed that the patch technique is only an improvement when oil is used as vehicle, since the pads soaked with aqueous solutions dried during over night exposure. Therefore, preparations in **oil** are applied at **20 µL plus pad**, whereas solutions in **H₂O** are applied at **50 µL without pad**.

- Solutions in H₂O: Apply **50µL** atop the EpiDerm™ tissue and gently spread.
- Solutions in oil: Soak the Finn pad by pipetting **20 µL** atop the pad and put it on the tissue.
- Suspensions in oil and other formulations: Soak "bottom" of Finn pad (tissue contact site) with 20 µL, turn the pad around and put it on the tissue.

5.5 Experimental Procedure

Day before testing

Upon receipt of EpiDerm™ kits (Europe: Tuesday afternoon), place assay medium and sealed 24-well plates containing tissues on agarose into refrigerator (4°C ± 2°C). Place the vial containing the MTT concentrate in the freezer (-20°C ± 5°C).

First day of testing

Note: It is essential that, before the test is started, tissues are incubated for **at least 1 hr** in assay medium into which they can release metabolites and debris accumulated during the shipment. This medium has to be replaced before the assay is started. For this important incubation 6-well plates have to be used with 0.9 mL medium per well. Alternative techniques (e.g. use of 24 well plates or incubation of all insets in a petri dish) had to be disapproved during this study).

1. Prewarm assay medium in a 37°C waterbath
2. Pipet 0.9 ml of assay medium into each well of sterile 6-well plates
3. Using sterile techniques transfer the inserts (be sure to remove all transport agar) into 6-well plates containing prewarmed assay medium. Any air bubbles trapped underneath the inserts should be released. Incubate for a minimum of 1 hr at 37°C, 5% CO₂ *Record incubation time in the MDS*
4. While tissues are in the incubator, for each test chemical, prepare a series of five concentrations according to **5.4.1**. *Record preparations in the MDS*
5. After (at least) 1 hr incubation transfer inserts into new 6-well plates prepared with new 0.9 mL assay medium per well and prewarmed in the incubator.
6. Per each test chemical use one EpiDerm™ kit (24 tissues): Twelve tissues are used in the **(-UVA) cytotoxicity** part and 12 in the **(+UVA) phototoxicity** part of

the test. Both parts of the test are dosed identically: Apply vehicle control (VC) plus 5 concentrations of the test chemical (C1 -C5), each on tissue duplicates according to **5.4.2**. Mark lids and plates to prevent from transposition errors.

+ UVA Plate 1

VC	C 1	C 2
VC	C 1	C 2

- UVA Plate 1

VC	C 1	C 2
VC	C 1	C 2

+ UVA Plate 2

C 3	C 4	C 5
C 3	C 4	C 5

- UVA Plate 2

C 3	C 4	C 5
C 3	C 4	C 5

Note: Covering the plate with self adhesive plastic film to prevent from toxic evaporation into neighbour wells is, in general, not necessary in this test. To be on the safe side, position the lowest concentration (C1) beside the vehicle control (VC).

- Once all tissues have been dosed, cover the plates with the lids and incubate over night (**18-24 hrs**) at 37°C, 5% CO₂.

Second day of testing

1. Remove 6-well plates from the incubator. Remove application pads. Irradiate (+UVA)-plates (covered with lids) for **60 min** with **1.7 mW/cm² (= 6 J/cm²)** at room temperature. Ventilate with fan to prevent condensation under the lid. Place (-UVA)-plates in the dark at room temperature.
2. While tissues are irradiated, prepare appropriate amount of new 6-well plates with 0.9 mL of fresh assay medium per well and prewarm in the incubator.
3. After UVA irradiation is completed, use wash bottle with sterile PBS and rinse each insert of the (+UVA) plates and (-UVA) plates. Then transfer all inserts to the new plates prepared in 2.
4. Incubate (+UVA)- and (-UVA) plates over night (18-24 hrs) at 37°C, 5 % CO₂

Third day of testing

1. Prepare MTT medium: thaw MTT concentrate in a water bath and dilute with MTT diluent. Spin down (300×g, 5 min) to remove any precipitate. Prewarm MTT medium (water bath) to 37°C.
2. For *each* test material prepare *one* 24-well plate with 300 µL prewarmed MTT medium per well. Label plates (lid and bottom) and transfer tissue inserts according to the plate design given below. Any air bubbles trapped underneath the inserts should be released.

VC	C 1	C 2	C 3	C 4	C 5	+UVA
VC	C 1	C 2	C 3	C 4	C 5	+UVA
VC	C 1	C 2	C 3	C 4	C 5	-UVA
VC	C 1	C 2	C 3	C 4	C 5	-UVA

4. Incubate 24 well plate 3 hours (37°C, 5 % CO₂). *Record start and stop time for MTT incubation in the MDS.*

Note: Deviations from 3 hour time for MTT incubation will result in different MTT readings. For consistency it is recommended that 3 hour MTT incubation time be adhered very strictly.

5. After incubation aspirate MTT medium (gently using a suction pump), refill wells with PBS and aspirate PBS. Repeat the procedure twice and make sure tissues are dry after the last aspiration. Transfer inserts to new 24 well plates.
6. For formazan extraction immerse the inserts by gently pipetting 2 mL extractant solution (isopropanol) **into** each insert. The level will rise above the upper edge of the insert, thus completely covering the tissue from both sides.
7. Seal the 24 well plate (e.g. with a zip bag) to inhibit isopropanol evaporation. *Record start time of extraction in the MDS.* Extract 2 hrs with shaking (~120 rpm) at room temperature.

8. After formazan extraction period is complete pierce the inserts with an injection needle (~ gauge 20 / 0.9 mm Ø) and allow the extract to run into the well from which the insert was taken. Afterwards the insert can be discarded. Place the 24-well plates on a shaker for 15 minutes until solution is homogeneous in colour.
9. Per each tissue transfer $3 \times 200\mu\text{L}$ aliquots* of the blue formazan solution into a 96-well flat bottom microtiter plate. For the 96 well plate, use **exactly the plate design given below** as this configuration is used in the EXCEL data spreadsheet. Read OD in a plate spectrophotometer at **570 nm, without reference filter**.# Alternatively, ODs can be read at **540 nm**.

* **Note:** In contrast to normal photometers, in plate readers pipetting errors influence the OD. Therefore, 3 formazan aliquots shall be taken from each tissue extract. In the data sheet these 3 aliquots will be automatically reduced to one value by calculating the mean of the three aliquots. Thus, for calculations from each single tissue only one single mean OD-value is used.

Note: Readings are performed without reference filter, since the "classical" reference filter often used in the MTT test (630 nm) is still within the absorption curve of formazan. Since filters may have a \pm tolerance in some cases the reference filter reduces the dynamics of the signal (OD) up to 40%.

Fixed 96 well-plate design (for OD reading in plate photometer, 3 aliquots per tissue)

VC	VC	C1	C1	C2	C2	C3	C3	C4	C4	C5	C5	+UVA
VC	VC	C1	C1	C2	C2	C3	C3	C4	C4	C5	C5	
VC	VC	C1	C1	C2	C2	C3	C3	C4	C4	C5	C5	
VC	VC	C1	C1	C2	C2	C3	C3	C4	C4	C5	C5	-UVA
VC	VC	C1	C1	C2	C2	C3	C3	C4	C4	C5	C5	
VC	VC	C1	C1	C2	C2	C3	C3	C4	C4	C5	C5	
tissue	tissue	tissue	tissue	tissue	tissue	tissue	tissue	tissue	tissue	tissue	tissue	
1	2	1	2	1	2	1	2	1	2	1	2	

5.6 Documentation

5.6.1 Method Documentation Sheet, MDS

The MDS allows to check the correct set up, calibration and function of the equipment as well as correct weights, applications etc. The MDS is designed as a paper document "in the spirit of GLP". Per each kit, make a hardcopy of the MDS, fill in and sign the requested information, starting the day prior to testing and ending after the test has been conducted.

final version 5 Novemb. 97	Standard Operation Procedure EpiDerm™ Phototoxicity Assay (model: Epi-200)	page 16 of 28
-------------------------------	---	------------------

Note (1): If several tests are performed per week, pipette verification (weighing H₂O on a balance) is only necessary once at the beginning of each week. Nevertheless, if **adjustable** pipettes are used the correct adjustment shall be checked and recorded in the MDS before each test.

5.6.2 Data Spreadsheet

The MS EXCEL 5.0 spreadsheet P-SPREAD.XLS is provided by ZEBET. Data files of optical densities (ODs) generated by the microplate reader are copied from the reader software to the Windows Clipboard and then pasted into the first map of the EXCEL spreadsheet in the fixed 96-well format given above (**Note:** Only 72 wells of the 96 wells are used!).

P-SPREAD.XLS consists of **two** maps, IMPORT, and SPREAD. The first map (*Import*) is used for pasting OD values (cursor position: **A20!**). The second map (*Spread*) does the calculations and provides a column graph of the results. In addition, entry all information requested (tissue lot-no., test material codes, date...) into this map.

In Phase III of the prevalidation study test chemicals are coded by BIBRA with a four digit code. To allow an easy allocation of the XLS files for statistical analysis after codes are broken, use the following file names **PGXXXX-Y.XLS** or **BDXXXX-Y.XLS** or **ZEXXXXX-Y.XLS** where XXXX stands for the 4-digit code number and Y stands for the number of the test run.

6 PREDICTION MODEL

The rules used to transform quantitative or qualitative data of a toxicological test into a prediction of a toxic potential or potency are called prediction model.

The prediction model is based on analysis of historical data of the maximum possible difference in the viability of identically treated EpiDerm™ tissues according to which any difference exceeding 30% has to be regarded significant ($p < 0.001$). Since the UVA irradiation has no cytotoxic effect itself, a phototoxic activity can be predicted if viability of tissues treated with identical test chemical concentrations differs by more than 30% in the irradiated and the non irradiated part of the test.

For each concentration of a test chemical, the mean OD of the tissue couple treated with this concentration is determined and expressed as relative percentage viability of the untreated vehicle controls. Identical calculations are performed for the (+UVA) part of the test and the (-UV) part of the test.

A chemical is predicted to have a phototoxic potential if one or more test concentrations of the (+UVA) part of the experiment reveal a decrease in viability **exceeding 30%** when compared with identical concentrations of the (-UVA) part of the experiment.

Prediction of phototoxicity is supported if, in addition, the (+UVA) induced reduction in tissue viability shows a dose response relationship.

Note: *It is a quite common observation for certain phototoxins, that a severe effect may be reduced again at higher doses. This is due to the UV absorbing properties of the chemical by which they act as UV filters if excessive doses are applied which remain on the stratum corneum.*


final version 5 Novemb. 97	Standard Operation Procedure EpiDerm™ Phototoxicity Assay (model: Epi-200)	page 17 of 28
-------------------------------	---	------------------

7 REFERENCES


1. Mosmann, T. (1983)
Rapid colorimetric assay for cellular growth and survival: Application to proliferation and cytotoxicity assays
J. Immun. Meth. **65**: 55-63.
2. Klausner, M., Kubilus, J., Ricker, H.A. and P.J. Neal (1995)
UVB irradiation of an organotypic skin model, EpiDerm™, results in significant release of cytokines.
The Toxicologist, *15* (1), *Soc. of Toxicology*.
3. Perkins, M.A., Osborne, R. and G. R. Johnson (1996)
Development of an *in vitro* Method for Skin Corrosion Testing.
Fundamental and Applied Toxicology **31**: 9-18.
4. Edwards, S. M., Donally T.A., Sayre, R.M. Rheins, L.A., Spielmann, H. and Liebsch. M. (1994):
Quantitative *in vitro* assessment of phototoxicity using a human skin model; Skin².
Photodermatol. Photoimmunol. Photomed **10**: 111-117
5. Rouget, R., Cohen, C. and A. Rougier (1994):
A reconstituted human Epidermis to assess cutaneous irritation, photoirritation and photoprotection *in vitro*. **In**: Alternative Methods in Toxicology, Vol. 10: *In vitro* Skin Toxicology - Irritation, Phototoxicity, Sensitization. **Eds.** A. Rougier, A. Goldberg, H. Maibach;
Mary Ann Liebert Publ., New York; pp. 141 - 149
6. Spielmann, H., Liebsch, M., Pape, W.J.W., Balls, M., Dupius, J., Klecak, G., Lovell, W.W., Maurer, De Silva, O., Steiling, W. (1995):
The EEC COLIPA *in vitro* photoirritancy program: results of the first stage of validation. **In**: Irritant Dermatitis: New clinical and experimental aspects. **Eds.** P. Elsner and H.I. Maibach
T. Karger Publ., Basel; pp. 256-264
7. Liebsch, M., Spielmann, H., Balls, M., Brand, M., Döring, B., Dupuis, J., Holzhütter, H.G., Klecak, G., L'Eplattenier, H., Lovell, W.W., Maurer, T., Moldenhauer, F., Moore, L., Pape, W.J.W., Pfannenbecker, U., Potthast, J., De Silva, O., Steiling, W., Willshaw, A. (1994):
First results of the EC/COLIPA validation project "*in vitro* phototoxicity testing". **In**: Alternative Methods in Toxicology, Vol. 10: *In vitro* Skin Toxicology - Irritation, Phototoxicity, Sensitization. **Eds.** A. Rougier, A. Goldberg, H. Maibach;
Mary Ann Liebert Publ., New York; pp. 243 - 254
8. Liebsch, M., Döring, B., Donnelly, T.A., Logemann, P., Rheins, L.A. and H. Spielmann (1995):
Application of the human dermal model Skin² ZK 1350 to phototoxicity and skin corrosivity testing.
Toxic. in Vitro Vol. 9, No. 4, 557 - 562

8 ANNEX A: Procedure Steps


First day of testing (Wednesday)



Second day of testing



Third day of testing (Friday)



final version 5 Novemb. 97	Standard Operation Procedure	page 19 of 28
	EpiDerm™ Phototoxicity Assay (model: Epi-200)	

9 ANNEX B: Methods Documentation Sheet

ASSAY No:.....

DATE:.....

XLS file name:.....

Test Chemical:.....

Kit receipt

EpiDerm kit received (day/date):	Day used:
EpiDerm Lot no.:	Production date:
Epi-100 Assay medium Lot no.:	Expiration date:
MTT concentrate Lot no.:	Date:
MTT diluent Lot no.:	Date:
MTT extractant Lot no.:	Date:
Booked in by (ID):	

PBS preparation (to be filled in only if PBS is prepared from concentrate / powder)

DPBS Lot no.:	Expiration date:	
Vol 10x DPBS:	Vol water:	Initial pH:
NaOH used to adjust pH:		Final pH:
HCl used to adjust pH:		Final pH:
Prepared by (ID):		

Incubator verification

Incubator #	CO ₂ (%)	Temperature (°C)	Check water in reservoir (✓)
			ID / date:

Pipette verification (triplicate weightings)

Note: Perform pipette verification only once per week and refer to it in all assays of this week. But: If adjustable pipettes are used, check correct adjustment daily and mark with (✓).

verification	0.9 mL	300 µL	200 µL	20 µL	50 µL
	H ₂ O weight (mg)				
1.					
2.					
3.					

ID / date:

final version 5 Novemb. 97	Standard Operation Procedure	page 20 of 28
	EpiDerm™ Phototoxicity Assay (model: Epi-200)	

Preparation of test chemical and application

		concentration (%)
stock preparation:		
1. dilution:		
2. dilution:		
3. dilution:		
4. dilution:		
5. dilution:		

vehicle used (H₂O / oil):

homogenisation technique:

solution / suspension:

application volume (µL):

applied + pad or - pad:

Time protocols

Procedure	Start	Stop
1 hr pre-incubation of tissues		
over night chemical application (incubator)		
3 hrs MTT incubation		
Formazan extraction		

ID / Date:

Check plate photometer filter (✓)


reading filter: 570 nm	
reading filter: 540 nm	

ID / Date:

Remarks



10 ANNEX C: EXCEL Spreadsheet

	chemical:			Bergamotteoil			solvent:		oil	
	tissue-lot no.:			1230 C			date of start:		07.05.97	
	exp. no.:			27			Irradiation:		6 J/cm ²	
							application time:		21 hours	
				mean						mean
control	0,920	0,902	0,887	0,903		control	1,117	1,168	1,185	1,157
	1,058	1,052	1,030	1,047			0,884	0,884	0,896	0,888
0,1	1,045	1,041	1,038	1,041		0,1	1,079	1,076	1,093	1,083
	0,982	0,951	0,943	0,959			1,060	1,054	1,088	1,067
0,316	1,051	1,055	1,012	1,039		0,316	0,975	0,970	0,979	0,975
	0,936	0,902	0,907	0,915			0,949	0,930	0,949	0,943
1	0,972	0,961	0,948	0,960		1	0,594	0,595	0,600	0,596
	0,867	0,861	0,854	0,861			0,622	0,625	0,625	0,624
3,16	0,969	0,969	0,958	0,965		3,16	0,312	0,310	0,313	0,312
	1,063	1,069	1,051	1,061			0,271	0,268	0,271	0,270
10	0,936	0,917	0,904	0,919		10	0,342	0,325	0,321	0,329
	0,917	0,933	0,912	0,921			0,281	0,280	0,280	0,280
	w/out UVA						with UVA			
	mean	Δ tissue [%]	% untreated control				mean	Δ tissue [%]	% untreated control	
control	0,975	14,74	100			control	1,022	26,28	100	
0,1	1,000	8,27	103			0,1	1,075	1,43	105	
0,316	0,977	12,72	100			0,316	0,959	3,34	94	
1	0,911	10,95	93			1	0,610	4,53	60	
3,16	1,013	9,44	104			3,16	0,291	14,33	28	
10	0,920	0,18	94			10	0,305	16,07	30	





final version 5 Novemb. 97	Standard Operation Procedure EpiDerm™ Phototoxicity Assay (model: Epi-200)	page 22 of 28
---	---	--------------------------------






11 ANNEX D: Positive Reference Data



12 ANNEX E: Negative Reference Data



13 ANNEX F: EPI-200 UVA-Sensitivity



final version 5 Novemb. 97	Standard Operation Procedure EpiDerm™ Phototoxicity Assay (model: Epi-200)	page 26 of 28
---	---	--------------------------------


The figures show four independent UV sensitivity experiments performed according to **5.3.1**. The dose of 6 J/cm² used in the EpiDerm™ Phototoxicity Test is not cytotoxic in any of the experiments. In addition, the dose of 6 J/cm² is comparable to doses used in animal tests and has proved to be sufficient to activate phototoxins.

final version 5 Novemb. 97	Standard Operation Procedure EpiDerm™ Phototoxicity Assay (model: Epi-200)	page 27 of 28
-------------------------------	---	------------------

14 ANNEX G Irradiance spectrum of the sun simulator

Note: The irradiance spectra of the SOL 3 and SOL 500 are nearly identical up to a wavelength of ~ 550 nm. In the longer wavelength range of visible light (> 550 nm - 700 nm) the SOL 500 irradiance decreases, whereas the irradiance of the SOL 3 remains at the same level. In the EU/COLIPA validation study this difference proved to irrelevant.

(spectrum kindly provided by Beiersdorf AG)



ADDENDUM (13 November 1997)

The spectrum (provided by Beiersdorf AG) shows the FORMAZAN absorption. It explains why the reference filter of 630 nm has been omitted. If the filter is not precise (e.g. 620 nm) the dynamics of the reading will be reduced by ~40% !!!

